

# Alaska Scientific Crime Detection Laboratory

---

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Archived

STATE OF ALASKA

DEPARTMENT OF PUBLIC SAFETY

SCIENTIFIC CRIME DETECTION LABORATORY

## BREATH ALCOHOL PROCEDURE MANUAL DataMaster DMT



---

Nita Bolz – Supervisor, Scientific Director

---

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### TABLE OF CONTENTS

#### SECTION 1 INTRODUCTION

- 1.1 Overview of the Alaska Breath Alcohol Program
- 1.2 Evidential Breath Test Instrument

#### SECTION 2 ADMINISTRATIVE PROCEDURES

- 2.1 Verification of Calibration Reports
- 2.2 DM Host
- 2.3 Laboratory Information Management System (LIMS)
- 2.4 Legal

#### SECTION 3 TECHNICAL PROCEDURES

- 3.1 Instrument Check-In
- 3.2 Instrument Calibration
- 3.3 Instrument Certification
- 3.4 Instrument Repair
- 3.5 Preventative Maintenance

#### SECTION 4 QUALITY ASSURANCE

- 4.1 External Dry Gas Ethanol Standards
- 4.2 Wet Bath Ethanol Simulator Solutions
- 4.3 Wet Bath Simulators
- 4.4 Barometer

#### SECTION 5 TRAINING PROVIDED

- 5.1 Breath Test Supervisor Training

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### 5.2 Breath Test Operator Training

#### **TABLE OF CONTENTS (continued)**

APPENDIX I	APPROVED ASSOCIATED EQUIPMENT
APPENDIX II	JUSTICE TRAX WORKING INSTRUCTIONS
APPENDIX III	BREATH ALCOHOL PROGRAM FORMS

# Alaska Scientific Crime Detection Laboratory

---

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### SECTION 1

#### INTRODUCTION

##### 1.1 OVERVIEW OF THE ALASKA BREATH ALCOHOL PROGRAM

The Scientific Crime Detection Laboratory (crime lab) is part of the Department of Public Safety and per 13 AAC 63.010 the commissioner will designate an employee of the crime lab as the scientific director of the blood and breath alcohol testing program. The scientific director is responsible for all aspects of the breath alcohol program including calibration and certification of breath test instruments; training and certification of breath test operators and breath test supervisors; and maintaining the scientific integrity of the breath test program.

The breath alcohol section staff falls under the supervision of the scientific director and as his/her designee(s) carry out the tasks associated with running the statewide breath testing program.

##### 1.2 EVIDENTIAL BREATH TEST INSTRUMENT

Per 13 AAC 63.020 and 13 AAC 63.030 the scientific director must approve a type of breath test instrument for use in the state of Alaska as well as maintain a list of all approved breath test instruments and associated equipment. The current breath test instruments approved for use in the state of Alaska are the DataMaster cdm and the DataMaster DMT both manufactured by National Patent Analytical Systems, Inc. The contents of this manual are specific to the DataMaster DMT. The approved associated equipment list is maintained in Appendix I of this manual.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

## SECTION 2

### ADMINISTRATIVE PROCEDURES

#### 2.1 VERIFICATION OF CALIBRATION REPORTS

Per 13 AAC 63.100 the accuracy of the calibration of the breath test instruments must be verified at least every 60 days. The documentation of this accuracy check is the verification of calibration report. The accuracy check performed for the verification of calibration report consists of a diagnostic check and five tests of the external dry gas ethanol standard. The verification of calibration procedure is incorporated into the instrument software under the direction of the scientific director. A verification of calibration is initiated by the instrument software in intervals of less than 60 days and whenever the external dry gas ethanol standard is changed. Breath test operators and breath test supervisors can also initiate a verification of calibration.

The instrument software is scheduled to initiate a verification of calibration on the following dates at 12:00 pm:

January 3<sup>rd</sup>  
February 21<sup>st</sup>  
April 11<sup>th</sup>  
May 30<sup>th</sup>  
July 18<sup>th</sup>  
September 5<sup>th</sup>  
October 24<sup>th</sup>  
December 12<sup>th</sup>

At the completion of a successful verification of calibration the instrument stores a copy of the written report to memory. If a verification of calibration does not complete successfully then the instrument will not allow a subject test to be initiated until a valid verification of calibration has been completed.

The instrument memory is periodically uploaded to a crime lab computer via a modem line to the DM Host software. The completed verification of calibration report is printed from the uploaded file and technically reviewed by the members of the breath test section prior to being signed and notarized by the scientific director. The signed copies are scanned and a copy is placed on the crime lab website and in the Laboratory Information Management System (LIMS). The original signed copy is retained in a file at the crime lab by the scientific director.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The dates listed above are not the only acceptable dates for verification of calibrations. The only requirement is that a verification of calibration be performed every 60 days when an instrument is in service.

### Technical Review of Verification of Calibration Reports

A technical review of each verification of calibration report is performed by a qualified member of the breath alcohol section prior to the report being reviewed and signed by the scientific director. The technical review includes the following components:

- Verification that the external dry gas ethanol standard cylinder lot number and expiration date is included in the list of approved cylinders kept at the crime lab.
- Verification that the five external standard results fall within  $\pm 0.005$  of the target value adjusted for barometric pressure.
- Verification that the standard deviation of the five external standard results is equal to or less than 0.0030.
- Ensure that each component of the diagnostic check reads "passed".
- Ensure no status messages are present.

The technical reviewer's initials and the date the review was performed are placed in the lower right-hand corner of the verification of calibration report indicating that the technical review was performed.

The scientific director completes a review that includes all components of a technical review and then signs the approved report in front of a notary. A verification of calibration is not a completed approved report until signed by the scientific director and notarized. The completed verification of calibration report certifies the instrument for continued evidentiary use in the State of Alaska.

## **2.2 DM HOST**

The DM Host software is provided by National Patent Analytical Systems, Inc. The DM Host software communicates, through modem lines, with the breath test instruments to upload the data written to memory. The DM Host software also enables the breath analysts to remotely access the instruments in the field for troubleshooting purposes. Breath analysts can view the technician screen remotely for troubleshooting purposes; however, the voltages and settings can only be adjusted when an instrument is at the crime lab.

### Uploading Data From Breath Test Instruments

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The DataMaster DMT writes subject tests, diagnostic tests, verifications of calibration, supervisor tests, linearity tests, status messages, calibration records, tank change records and software update records to the instrument memory. The DM Host software enables the breath test section to manually upload the data collected by the breath test instruments in the field. The data is then stored in the DM Host database which is housed on a secure server at the crime lab. Once the data is uploaded to the crime lab it is no longer available in the DataMaster DMT memory.

All data retrieved by DM Host from the breath test instruments is considered an official record and that data is not altered by the breath alcohol section staff. The data collected can be searched, reports generated and printed. Breath test instruments in the field will be uploaded periodically.

### 2.3 LABORATORY INFORMATION MANAGEMENT SYSTEM (LIMS)

The LIMS currently in use at the laboratory is Justice Trax. The Justice Trax system is used by the breath alcohol program to electronically file the paperwork associated with running the breath test program.

The Justice Trax program is used to track breath test operators and breath test supervisors along with their training records. The program also tracks the status and location of each instrument along with any paperwork associated with that instrument.

The working instructions for the operation of Justice Trax are located in Appendix II of this manual.

### 2.4 LEGAL

One of the duties of the breath alcohol program is to work with the legal community to provide testimony and information about the breath alcohol program, the evidentiary breath testing instrument and alcohol impairment.

#### Discovery

Discovery requests are handled by the scientific director or his/her designee. The breath alcohol program follows the laboratory discovery policy. For more information on the laboratory discovery policy refer to the Laboratory Policy Number 10: Disclosure of Examination Results.

#### Testimony

# Alaska Scientific Crime Detection Laboratory

---

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The breath alcohol analysts provide testimony on the breath alcohol program, the evidentiary breath testing instrument and alcohol impairment. This testimony is provided in state, municipal and federal courts upon request and approval by the section supervisor.

### Opinions Issued

On occasion the breath alcohol program is asked to provide written opinions on various aspects of breath alcohol or impairment related topics. These opinions are generated as memos either from the breath alcohol section as a whole or from the individual breath alcohol analyst depending on the specific topics covered. Opinions are not issued as reports and therefore they are not technically reviewed. Opinion memos are scanned into the LIMS after completion.

### Department of Motor Vehicles

The Department of Motor Vehicles (DMV) holds administrative proceedings regarding impaired driving cases and can request information regarding these cases. Requests from the DMV are handled following the laboratory discovery policy. Occasionally, expert testimony is provided for these hearings upon request and approval by the section supervisor.



# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### SECTION 3

#### TECHNICAL PROCEDURES

##### 3.1 INSTRUMENT CHECK-IN

When an instrument is initially received by the laboratory from the manufacturer a DataMaster Acceptance form is filled out by the receiving analyst. This information is then stored in the LIMS.

When instruments arrive at the laboratory from the field, the breath test instrument will follow a check-in procedure. The check-in procedure includes updating the instrument's status in the LIMS, initiating a verification of calibration, printing the calibration factors and uploading the instrument to retrieve all information stored in memory.

A qualified breath analyst will then perform a visual inspection and evaluate the instrument for repair. If a repair is required and can be performed at the laboratory, the analyst will complete the repair process. If the repair requires factory assistance, the instrument will be sent to National Patent Analytical Systems, Inc. for repair. The check-in process along with any repair(s) performed should be documented on the check-in form which is scanned into the LIMS.

The breath test supervisors communicate with the breath alcohol section about the status of their agency's instrument utilizing a change in instrument status form. This form includes information on when an instrument was placed in service, taken out of service and what problems, if any, reported with the instrument.

All correspondence from breath test supervisors recorded on the change in instrument status forms will be documented in the LIMS by the breath alcohol section.

##### 3.2 INSTRUMENT CALIBRATION

Per 13 AAC 63.100 the scientific director has designated the DataMaster DMT can be calibrated by either National Patent Analytical Systems, Inc. or a qualified breath analyst at the crime laboratory. The DataMaster DMT calibration is a single point calibration. Once instruments are calibrated they are not recalibrated unless previously stored calibration factors are no longer accessible by the DataMaster DMT rendering the instrument inoperable or a breath analyst determines the instrument is not functioning properly.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### Calibration Procedure

The crime lab performs calibration of the DataMaster DMT using a wet bath ethanol simulator solution and a blank simulator solution. The blank simulator solution is 500 mL of deionized water.

Wet bath ethanol simulator solutions used by the crime lab are commercially prepared, traceable to National Institutes of Standards and Technologies (NIST), and come accompanied by a certificate of analysis. These solutions have a concentration that yields an ethanol vapor concentration at 34°C of 0.10 g/210L on a breath test instrument. The certificates of analysis for simulator solutions are scanned into the LIMS.

The crime lab uses simulators manufactured by Guth Laboratories that are NIST traceable. They are returned to Guth Laboratories annually for calibration and maintenance. Certificates of calibration for the simulators are scanned into the LIMS.

The DataMaster DMT software allows the breath analyst to select the number of replicate samples of the simulator that the instrument requires for calibration. This option is found under the DMT Icon Setup menu for calibration. This should be set to run three replicates for both the blank and 0.10 g/210L simulator solutions.

When a qualified breath analyst determines that an instrument requires recalibration the analyst should create a calibration request in the LIMS and assign the request. A new bottle of simulator solution should be poured into a simulator and the date that the bottle was opened indicated on the label of the simulator. A new bottle of simulator solution should be used each day that instruments are calibrated. The analyst should ensure the instrument and simulators have been turned on and warmed up prior to performing calibration. The lot number and expiration date of the simulator solution and serial number and calibration date of the simulator used is documented in the LIMS.

The DataMaster DMT instrument software prompts the breath analyst to attach both the blank and 0.10 g/210L ethanol simulator solutions to the breath instrument's calibration port at the appropriate times during the calibration sequence. The breath analyst should follow the instrument prompts. When the calibration process has been completed the breath analyst should print a copy of the calibration factors from the new calibration. The calibration factors should be scanned into the laboratory calibration request in the LIMS. Each laboratory calibration process is technically reviewed by another qualified analyst. When the technical review has been completed a certificate of calibration will be generated.

When calibration is performed by the manufacturer it is returned to the laboratory with a certificate of calibration. This certificate should be scanned into the LIMS under a factory calibration request. Every time the instrument is calibrated (whether by the

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

factory or crime lab) the instrument must have the calibration verified with a certification procedure.

### 3.3 INSTRUMENT CERTIFICATION

As per 13 AAC 63.100 every breath test instrument must have the calibration verified prior to being certified for use. This must be performed by the scientific director or a qualified person designated by the scientific director. Instrument certification is the process that the breath alcohol section uses to verify that the individual breath instrument's calibration has been verified. The scientific director has designated that instrument certifications can be performed by qualified breath analysts from the crime lab.

An instrument is required to have a certification process performed after the instrument is received from the manufacturer, after a repair that has an effect on the analytical components of the instrument, and as part of preventative maintenance. Certification is not necessary after software updates unless the update has an effect on the analytical components of the instrument.

The DataMaster DMT certification process includes voltage checks, linearity testing, functional tests, minimum volume determination and a non-drinking subject test.

#### Voltage Checks

The DataMaster DMT has digital potentiometers that can be viewed and adjusted by the breath analysts in the technician screen located under the DMT Icon. During certification the technician screen is checked by the analyst to ensure the voltages displayed are within normal operating tolerances published by the instrument manufacturer. The technician screen is also the location of the barometer setting, the radio frequency detection sensitivity and the breath / air volume and flow rate readout.

The barometer reading of the breath test instrument should be checked with the NIST traceable barometer located at the crime lab and the setting should be adjusted if the instrument barometer differs by more than +/- 0.10 inches of mercury. The barometer used by the crime lab for this comparison must be purchased from an ISO 17025 certified supplier, be NIST traceable and come accompanied with a certificate of calibration. The barometer should be sent out for recalibration or replaced when the manufacturer calibration expires. The certificates of calibration for the barometer will be scanned into the LIMS.

The radio frequency sensitivity is adjusted from the technician screen. Prior to adjusting the sensitivity the analyst should ensure that all radios are turned off in the vicinity of the

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

instrument. The initial adjustment should be performed by hitting the set button. This will trigger a series of beeps as the instrument sets a default baseline level from which it measures radio frequency.

The volume adjustment ensures that the DataMaster DMT is appropriately reading breath volume. The breath analyst should use a graduated syringe to provide a 1.5 L breath sample to the DataMaster DMT while in the technician screen. The breath analyst should ensure that the volume reading in the technician screen reads within  $\pm 0.1$  L or the volume reading should be adjusted using the set screw inside the instrument.

### Linearity Testing

The linearity test checks the DataMaster DMT's ability to accurately read ethanol at a range of levels. This is accomplished by using the linearity test function. A linearity test should be conducted using external dry gas ethanol standards at the following levels: 52 ppm, 103.5 ppm, 208 ppm, 260.5 ppm, 521 ppm and 781 ppm. These levels are the equivalent of 0.020, 0.040, 0.080, 0.100, 0.200, and 0.300 g/210L respectively as measured by a breath test instrument. Each external dry gas ethanol standard is analyzed five times. The allowable range for the external dry gas ethanol standard is  $\pm 0.005$  or 3% , whichever is greater, of the target value adjusted for barometric pressure. The standard deviation of the five values must not be greater than 0.0030. The linearity test in the DataMaster also gives an R-squared value and the requirement for this is at least 0.999.

### Functional Tests

The DataMaster DMT generates status messages to inform the operator when certain situations occur. As part of the certification process some of these status messages are triggered to ensure the instrument is functioning properly. It is not possible to trigger every status message the DataMaster DMT can generate because some status messages require there to be a functional problem with the instrument. The following is a list of the status messages that are triggered during the certification process and their purpose.

- Invalid-An invalid status message is generated when the DataMaster DMT detects a negative slope in the subject's breath alcohol profile. This can be an indication of "mouth alcohol" or a high concentration of ethanol in the subject's mouth. This status message can be triggered by swishing ethanol containing mouthwash in the analyst's mouth prior to providing a breath sample to the instrument.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

- **Incomplete-** An incomplete status message is generated when the subject fails to provide a breath sample that meets the minimum requirements of the instrument within the two minutes allotted. To trigger this status message the instrument should be set up for a subject test and less than a minimum sample provided.
- **Interference Detected-** An interference detected status message is generated when the DataMaster detects the presence of a substance other than ethanol in the sample provided. Two of the more common interfering substances that could be found in a subject breath sample are acetone and isopropanol. To generate this status message an ethanol wet bath simulator solution containing 200 uL of acetone should be provided to the instrument through the breath hose. This is then repeated using an ethanol wet bath simulator solution with 200 uL of isopropanol. The solutions used for interference testing are made at the crime lab using a commercially prepared wet bath ethanol simulator solutions and adding 200 uL of either acetone or isopropanol. The lot number, ethanol level and expiration date of the simulator solution and the serial number and calibration date of the simulator used should be documented in the certification paperwork. The wet bath ethanol simulator solutions used for interference detected testing should be replaced three months after preparation.
- **Radio Frequency Interference Detected-** A radio frequency interference detected status message is generated when the DataMaster DMT detects the presence of radio frequency in the vicinity of the instrument. This status message is generated by setting up a subject test sequence and keying a radio in the room with the DataMaster. The instrument should generate an RFI status message when a radio is keyed in the room with the DataMaster DMT.

### Minimum Volume Determination

Volume checks ensure that the instrument is requiring the appropriate minimum breath volume. The DataMaster DMT is set to accept breath samples at a minimum volume of approximately 1.5 L.

Volume checks are performed using a graduated syringe. The breath analyst should provide samples in increments of 0.1 L through the breath hose to determine the minimum volume that the instrument being certified will accept. The requirement for an

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

instrument to be placed in service is that the minimum volume be greater than or equal to 1.0 L and less than or equal to 1.7 L.

### Non-Drinking Subject Test

The final step in an instrument certification process is a non-drinking subject test. The breath analyst should insert a 208 ppm ethanol (0.080 g/210L) external dry gas ethanol standard cylinder and follow the tank change procedure. This will initiate a verification of calibration. The analyst should then initiate a subject test using default software options. The subject test printout should be attached to the certification paperwork and the verification of calibration should be technically reviewed and signed by the scientific director. The verification of calibration report is not part of the certification paperwork.

If an instrument fails any portion of the certification process the instrument should be evaluated for repair or recalibration by the breath analyst. The certification form filled out for the initial certification paperwork should be scanned into the calibration request image module if recalibration was required. If a repair is required the initial certification paperwork should be scanned into LIMS where appropriate.

Once an instrument has completed the certification process the breath analyst should assign the certification request in the LIMS. The certification paperwork should be scanned into the request and the request milestone rolled to draft complete. The instrument certification must be technically reviewed by another qualified analyst prior to the instrument being sent out for use. This technical review is documented in the LIMS. Once the certification is technically reviewed a certificate is generated stating the instrument is certified for use.

### **3.4 INSTRUMENT REPAIR**

Once an instrument is certified for use and placed in the field it is the responsibility of the breath test supervisors to keep the crime lab informed of operational problems with the instruments. Breath test operators are trained to contact their breath test supervisor with problems who in turn can contact the crime lab if they are unable to correct the problem themselves.

The crime laboratory has several options for troubleshooting an instrument and the steps used in instrument troubleshooting are generally left up to the individual breath analyst.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

Breath analysts can work with breath test supervisors and/or remotely access the DataMaster DMT using the DM Host software. If the instrument is deemed to need repair by the crime lab, the breath test supervisor will take the instrument out of service and return it to the crime laboratory for repair.

Instruments arriving at the crime laboratory for repair should follow the check-in procedure and then a qualified breath analyst should assess the instrument to determine what repairs are needed, if any. This troubleshooting is performed using training provided by the instrument manufacturer and often involves direct contact with the manufacturer.

When the repair is completed, either by the factory or the crime lab, the instrument must undergo a certification prior to being placed in service at an agency. Exchanging the external accessories such as printers, keyboards and breath hoses do not require recertification because they do not alter the analytical components of the instrument. Software updates only require certification if they alter the analytical components of the instrument.

All repairs performed by the crime lab should be documented on the check in sheet. Information about repairs completed by National Patent Analytical Systems, Inc. is scanned into the LIMS as part of the check in paperwork.

### 3.5 PREVENTATIVE MAINTENANCE

In order to ensure that instruments in the field remain in good working order the crime lab has instituted a preventative maintenance schedule. This ensures that every instrument is returned to the crime lab for a visual inspection and recertification at least every four years. If an instrument is returned to the laboratory for repair the preventative maintenance will be performed at that time and the four year period will be restarted. Preventative maintenance consists of a check-in procedure; visual inspection and repairs as needed; and certification.



# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### SECTION 4

#### QUALITY ASSURANCE

##### 4.1 EXTERNAL DRY GAS ETHANOL STANDARDS

Dry gas ethanol standards are used as external standards with instruments in the field and for linearity testing as part of the certification process performed at the crime lab. External dry gas ethanol standards are purchased from an ISO 17025 certified manufacturer, are NIST traceable, and come with a certificate of analysis. An accuracy check will be performed on each new external dry gas ethanol standard lot received by analyzing it five times on a certified DataMaster. The results of this accuracy check must fall within  $\pm 0.005$  or 3%, whichever is greater, of the target value adjusted for barometric pressure and the standard deviation must not be greater than 0.0030. The certificate of analysis and the results of the accuracy check are then scanned into the LIMS.

##### 4.2 WET BATH ETHANOL SIMULATOR SOLUTIONS

Commercially prepared wet bath ethanol simulator solutions are used for calibration of breath test instruments and for interference detected functional checks. Ethanol simulator solutions are purchased from a NIST traceable supplier and come with a certificate of analysis. When a new lot of simulator solution is received an accuracy check will be performed on one bottle from each lot. The accuracy check will consist of analyzing one bottle five times on a certified DataMaster. The results should read within  $\pm 0.005$  or 3%, whichever is greater, of the expected value and the standard deviation must not be greater than 0.0030. The certificate of analysis and accuracy test will be scanned into the LIMS.

##### 4.3 SIMULATORS

The crime lab uses simulators manufactured by Guth Laboratories that are NIST traceable and come with a calibration certificate. Simulators are returned to Guth Laboratories annually for calibration and maintenance. Certificates of calibration for the simulators are scanned into the LIMS.

##### 4.4 BAROMETER

The barometer used for checking and adjusting the barometric pressure is purchased from an ISO 17025 certified manufacturer, NIST traceable and comes with a certificate of calibration. The



# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

barometer is sent out for calibration or replaced when the manufacturer calibration expires. Certificates of calibration are scanned into the LIMS.

### SECTION 5

#### TRAINING PROVIDED

##### 5.1 BREATH TEST SUPERVISOR TRAINING

One of the duties of the breath alcohol section is to train and certify breath test supervisors. Per 13 AAC 63.070 a breath test supervisor must complete a course approved by the scientific director. The course content is written and taught by the breath alcohol section and approved by the scientific director.

To be certified by the scientific director as a breath test supervisor the applicant must obtain a score of 75 percent or higher on a written examination covering the content of the course. The examination is written by the breath alcohol section and approved by the scientific director. Upon successful completion of the breath test supervisor course a list is provided to the Alaska Police Standard Council and a certificate is issued indicating the breath test supervisor is authorized to teach breath test operator courses.

The breath test supervisor certification expires at 11:59 p.m. on December 31 of the third year after issued. Per 13 AAC 63.080 to recertify as a breath test supervisor the applicant must complete a recertification course approved by the scientific director. This course is also taught by the breath alcohol section.

##### 5.2 BREATH TEST OPERATOR TRAINING

Breath test operators are trained by breath test supervisors in accordance with 13 AAC 63.050 and 13 AAC 63.060. The breath alcohol section provides the approved training program and examination to breath test supervisors for this training.

To certify a breath test operator the breath test supervisor should complete a Request for Certification of Breath Test Operators on the DataMaster DMT form to the crime lab indicating that the applicant has successfully completed the certification or recertification course. The breath alcohol section will then issue or renew an operator certificate and operator number indicating that the applicant is certified by the scientific director as a breath test operator.

Per 13 AAC 63.050 a breath test operator certificate expires three years from the date issued.

# Alaska Scientific Crime Detection Laboratory

---

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### APPENDIX I

#### APPROVED ASSOCIATED EQUIPMENT

##### External Dry Gas Ethanol Standards

External dry gas ethanol standards are obtained from the breath alcohol section. External dry gas ethanol standards used by the breath alcohol section are purchased from an ISO 17025 certified manufacturer, are NIST traceable, and come accompanied by a certificate of analysis. Each lot received by the crime lab is approved for use by the breath alcohol section after an accuracy test is performed on the lot. For this reason only external dry gas ethanol standards provided by the breath alcohol section may be used with the evidential breath testing instruments.

##### Mouthpieces

Mouthpieces are obtained from the breath alcohol section or the Department of Public Safety Supply. The current mouthpieces used are part number GM-1000-100 purchased from Guth.

##### DataMaster Parts and External Accessories

All parts used for instrument repair are purchased from National Patent Analytical Systems, Inc. and are approved for use with the DataMaster DMT.

External accessories used with the DataMaster DMT include the breath hose and USB keyboard. Replacement parts are obtained from the breath alcohol section.

##### External Printers

The first external printer is provided by the breath alcohol section. It is an HP LaserJet Professional model number P1606dn. Replacement printers, paper and toner are the responsibility of the individual agency. Any compatible toner cartridges may be used with the printer. Replacement printers must be approved by the breath alcohol section prior to purchase.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active


or use to ensure they are compatible with the DataMaster DMT. Any compatible external printer can be used with the DataMaster DMT.

## APPENDIX II

### JUSTICE TRAX WORKING INSTRUCTIONS

#### LOG IN

Your log in for the breath alcohol database is different from your log in for the laboratory database. You will know you are in the breath alcohol database because the splash screen has been removed. Your login name will be your last name.



#### CASE NUMBERS

The case mask for breath alcohol cases in the LIMS database all begin with the letter B. Cases numbers are named for instrument serial numbers. Instrument serial numbers follow the layout


# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

of 100XXX. Case numbers are in the format B-100XXX. Quality assurance documents are also kept as a case file. The case is B-QA2011. You can't use the case mask to open quality assurance cases. The case name will change each year. The format will still follow B-QAYYYYY. In the image module of the QA documents file will be four folders containing QA information for simulators, simulator solutions, dry gas ethanol standards (sorted by certification standards and standards for field use) and the barometer for the associated year. Opinions such as retrograde memos, etc., are also kept as a case file. The case name will change each year. The format will still follow B-OPINIONYYYY. You can't use the case mask to open opinion cases. The image module of the Opinion documents are kept in two folders: section opinions and analyst opinions.



### EVIDENCE

Each case will have an instrument assigned to it as evidence. The evidence will be titled the instrument serial number.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The screenshot shows a software window titled "Case B-100354" with a sub-header "100354". The window has a tabbed interface with tabs for "Case Info", "Agency", "Offense", "Individuals", "Evidence" (which is active), "Requests", and "Local Data". Below the tabs, a text field contains "100354 Datamaster DMT (NPAS John Fusco 01/19/2011)".

The main form area is divided into two columns. The left column has a header with radio buttons for "Genealogy" (selected), "Container", and "Storage Location", followed by "View" and "Edit" radio buttons. The left column contains fields for "Description" (Datamaster DMT), "Notes", "Evidence Type", "Source", "Intended Disp" (NPAS), "Other ID" (with a "Hold Evid at the lab" checkbox), and "Origin". The right column contains fields for "NCIC Evid Lbl", "Submit. Agency" (NPAS), "Badge Agc. Rep" (Fusco, John), "Kit", "Location" (Breath Lab), "Container", "Inherit", "Evidence No" (100354, with a barcode icon), and "Extraction type".

At the bottom left, there are status indicators: a red square for "Hold" and a green square for "Returned". At the bottom right, there is a "Submit for Analysis" checkbox and a "Hide Evidence Detail" checkbox. The bottom of the window features "OK", "Cancel", and "Apply" buttons.

The instrument will be related to each request and the location of the instrument will be tracked using the chain of custody. Each instrument will have a database barcode sticker on the side panel for use in evidence transport. There will also be a binder containing agency labels that can be used to track the chain of custody when an instrument is sent into the field.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

Case B-100354 100354

Case Info Agency Offense Individuals Evidence Requests Local Data

Req #

**New Request for Analysis**

**Requesting Party Information**

Agency: Breath Program Request Date: 02/18/2011

Badge Rep: Program, Breath

**Request Information**

Lab: Breath

Department: DataMaster DMT

Service: DataMaster Acceptance Analyst: O'Bryant, Colleen S

Due Date: Reason: Complexity:

**Notes**

Requester:

Assignor:

OK Cancel SOP

Request Status: (C)anceled (P)ending (R)eleased (IP) In Progress

OK Cancel Apply Hide Evidence Detail

When working with new requests you need to relate the evidence (instrument) to your request and update the chain of custody. The breath test section is an agency representative (people) while National Patent Analytical Systems, Inc. and other law enforcement agencies are storage locations.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

**Evidence Transfer**

From: National Patent Analytical Systems

Time: 02/18/2011 10:12 AM VIA [dropdown] Note [text area]

To: Breath Test Section,

Time: 02/18/2011 10:12 AM VIA [dropdown] Note [text area]

Then To: [text area]

Changes made to the Transfer information above will apply only to newly selected evidence.

**Evidence to Transfer** Scan the evidence to be included in this transfer.

Lab Agencies Ca	Container	Sub. #	Description	From	Time
-----------------	-----------	--------	-------------	------	------

☐ 1st Transfer Receipt ☐ 2nd Transfer Receipt ☐ Returned ☐ Chain of Custody Reports

Apply Clear Close

**Select Storage Location**

**Storage**

- Breath
  - Agencies
  - B-FedEx
  - Breath Lab
  - National Patent Analytical Systems

OK Cancel


# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011


Version: BAPM2011 R0  
Status: Active

### TYPES OF REQUESTS



Req #	Service	Status	Req Date	Due Date	Rel Date	Milestone	Assigned To
0006	DataMaster Acceptance	IP	02/14/2011	//	//	Findings Entered	O'Bryant, Colleen S
0007	Certification	IP	02/14/2011	//	//	Findings Entered	O'Bryant, Colleen S
0008	Check-in	IP	02/14/2011	//	//	Findings Entered	O'Bryant, Colleen S
0009	Preventative Maintenance	C	02/14/2011	02/13/2015	//	Unassigned	
0010	Certification	IP	02/14/2011	//	//	Findings Entered	Barnett, Brandi
0011	Factory Calibration	IP	02/14/2011	//	//	Findings Entered	NPAS
0012	Check-in	IP	02/14/2011	//	//	Findings Entered	O'Bryant, Colleen S
0014	Certification	IP	02/14/2011	//	//	Findings Entered	Barnett, Brandi
0015	Laboratory Calibration	IP	02/14/2011	//	//	Findings Entered	Barnett, Brandi

DataMaster Acceptance – This is the request used when the instrument is received after purchase. First add the Breath Section as an agency for the case. Under the agency tab right click and choose add agency. The Breath Program is the agency and the instrument number is the agency case number. Check the box making this the primary agency.



Case B-100362

100362

Case Info Agency Offense Individuals Evidence Requests Local Data

Related Agency Agency's Case No. County

Breath Program 100362

**New Related Agency for Agency's Case No. B-100362**

Agency Name: Breath Program ☒ Primary Agency

Agency's Case No.: 100362 ☐ Use ACN Mask

OK Cancel

Note: ☒ Indicates primary Agency.

OK Cancel Apply ☐ Hide Evidence Detail



# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

Create a new DataMaster Acceptance request. The requesting agency and agency representative is Breath Program. The lab is breath and the department is DataMaster DMT. Select DataMaster Acceptance as the service and yourself as the analyst. Press OK and this will create your request.

**New Request for Analysis**

**Requesting Party Information**

Agency: Breath Program Request Date: 02/23/2011

Badge Rep: Program, Breath

**Request Information**

Lab: Breath

Department: DataMaster DMT

Service: DataMaster Acceptance Analyst: Barnett, Brandi

Due Date: / / Reason: Complexity:

**Notes**

Requester:

Assignor:

OK Cancel SOP

Relate the instrument (evidence) to your request and update the chain of custody to show the instrument as being at the laboratory and accepted by the section.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The screenshot shows the 'Evidence Transfer' window. It has a 'From' field with 'National Patent Analytical Systems', a 'Time' field with '02/18/2011 10:12 AM', a 'VIA' dropdown, and a 'Note' field. The 'To' field is 'Breath Test Section,'. Below this is another 'Time' field with '02/18/2011 10:12 AM', a 'VIA' dropdown, and a 'Note' field. A 'Then To' field is at the bottom. A message bar states: 'Changes made to the Transfer information above will apply only to newly selected evidence.' Below this is a table titled 'Evidence to Transfer' with columns: Lab Agency's Ca, Container, Sub. #, Description, From, and Time. The table is currently empty. At the bottom are checkboxes for '1st Transfer Receipt', '2nd Transfer Receipt', 'Returned', and 'Chain of Custody Reports', along with 'Apply', 'Clear', and 'Close' buttons.

Use the DataMaster DMT Acceptance Form to check the instrument in. Once this form has been completed scan it into the image module and put it in the DataMaster Acceptance request folder. The calibration certification from the factory should also be scanned into the DataMaster Acceptance folder in the imaging module.

The screenshot shows two windows. On the left is the 'LIMS - Imaging' window with a file tree showing folders like 'Case Images [2]', 'Evidence', 'Individuals', and 'Requests'. Under 'Requests', there are several sub-items including '0006: DataMaster Acceptance [1]', '02/14/2011 - DATAMASTER ACCEPTANCE FORM\_0', '0007: Certification [1]', '0008: Check-in [1]', '0009: Preventative Maintenance [0]', '0010: Certification [1]', '0011: Factory Calibration [1]', '0012: Check-in [1]', '0014: Certification [0]', and '0015: Laboratory Calibration [1]'. On the right is the 'QSWNV102610869.pdf - Adobe Reader' window showing the 'DataMaster DMT Acceptance Form'. The form includes the Alaska Scientific Crime Detection Laboratory logo, the title 'DataMaster Acceptance Form - DMAF2011R0', and fields for 'Issuing Authority - Nita Bolz, Scientific Director', 'Approval date: 02/17/11', 'Effective date: 02/17/11', and 'Archive date: Active'. The form also has fields for 'Instrument # 100XXV', 'Date: 2/14/11', and 'Analyst: C. H. Lee'. A list of 8 tasks is provided, with checkboxes for 'Visual Inspection', 'Upload', 'Barcode Instrument', and 'Verification of Calibration - Date: 2/14/11'. The 'State of Alaska Property ID#' is listed as '10140909ABC'.


# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

Use the edit findings function; enter a value into findings for the evidence item by right clicking and adding a result. (A space in the results field will suffice as a result.) No analytical result is needed. After you have entered the findings, press ok, it will ask if you want to create your cascading request of “Certification”, press yes and close the edit findings window. This report is complete at draft complete. A report is not generated. After the report is rolled to draft complete a cascading request of “Certification” will be created.



Check – In – This is the request used when an instrument has been brought back to the laboratory from the field. If you are checking in an instrument, create the request and assign it to yourself. You should relate the evidence (instrument) to your request and update the chain of custody for the instrument to indicate where it has come from and that it is back at the laboratory. You should use the information from the change in instrument status form, if the returning agency sent one, to update the chain of custody and Check-In Form. Use the DataMaster Check-In Form to document the reason for repairs needed and the check in procedure. A verification of calibration will be performed during the check in procedure. It needs to be technically reviewed, signed by the scientific director and notarized. The VOC will be scanned into the image module and placed in the VOC folder. This folder will need to be created if it has not been already. You create a folder by right clicking on the case images file in the image module and created a subfolder. Name this folder VOC.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The screenshot shows two windows. The left window is 'LIMS - Imaging' with a tree view containing 'B-TEST', 'Case Images [2]', 'Evidence', 'Individuals', and 'Requests'. The right window is 'QSWNV0116C339.pdf - Adobe Reader' displaying a 'DataMaster Check In Form - DMCIF2011R0'. The form includes fields for 'Instrument # 100xxx', 'Date: 2/14/11', and 'Analyst: C. Brown'. It lists seven steps: 1. Reason for instrument return: 'Persistent Ambient Fail'; 2. Change instrument status in JusticeTrax Breath Alcohol Database (checked); 3. Visual Inspection (checked); 4. Verification of Calibration - Date: 2/14/11; 5. Print calibration factors (attached); 6. Upload (checked); 7. Troubleshooting/Repair: 'Unable to duplicate problem. Performed VOC and ID supervisor tests'.

### Factory Repair Required

If the instrument needs to be returned to the factory for repairs note this on the check in form. Then update the chain of custody to show the instrument has left the laboratory and is going to the factory. Do not roll the request to draft complete until after the instrument has returned from the factory.

The 'Evidence Transfer' dialog box shows transfer details from 'Breath Test Section' to 'National Patent Analytical Systems' on '02/24/2011 09:19 AM' via 'FedEx'. It includes a table for 'Evidence to Transfer' with columns: Lab Agency's Ca, Container, Sub. #, Description, From, and Time. At the bottom, there are checkboxes for '1st Transfer Receipt', '2nd Transfer Receipt', 'Returned', and 'Chain of Custody Reports', along with 'Apply', 'Clear', and 'Close' buttons.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

To ship the instrument print an evidence barcode label to go on the outside of the box for ease of evidence transfer. You can do this from the evidence tab by clicking on the barcode icon and choosing generate.

The screenshot shows the 'Case B-TEST' application window. The 'Evidence' tab is selected, displaying various fields for evidence tracking. A context menu is visible over the 'Evidence No' field, offering actions like 'Register', 'Generate', 'Delete-Unregister', and 'Generate Batch of Evidence Barcodes'. The 'Generate' option is highlighted, corresponding to the instruction in the manual to click the barcode icon and choose 'generate'.

After the instrument has come back from the factory (if factory repair needed), update the chain of custody to reflect it has returned. Then scan any paperwork for factory repairs into the check in request in the image module. If the factory calibrated the instrument, create a Factory Calibration request and scan the paperwork into the image module. The Factory Calibration request is described in detail below.

### After Factory Repair, or if none was needed


Using the edit findings function, enter a value into findings for the evidence item in the check in request. (A space in the results field will suffice as a result.) No analytical result is needed. This report is complete at draft complete. A report is not generated. After the report is rolled to draft complete a cascading request of "Certification" will be created. The same analyst that signed the instrument in and related the request to themselves does not have to be the analyst that rolls the request to draft complete.

# Alaska Scientific Crime Detection Laboratory


## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active



Factory Calibration - When the instrument is calibrated at the factory it will return with a calibration certificate. Create a factory calibration request and select yourself as the analyst. Relate the instrument to the request.



Scan the factory calibration certificate into the factory calibration folder in the image module. Then using the edit findings function, enter a value into findings for the evidence item by right clicking and adding a result. (A space in the results field will suffice as a result.) No analytical result is needed. This report is complete at draft complete. A report is not generated.

Laboratory Calibration – When an instrument is calibrated at the laboratory a Calibration Form is generated. Create a laboratory calibration request and select yourself as the analyst. Relate the instrument to the request.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

**New Request for Analysis**

**Requesting Party Information**

**Agency** Breath Program **Request Date** 02/24/2011

**Badge Rep** Program, Breath

**Request Information**

**Lab** Breath

**Department** DataMaster DMT

**Service** Laboratory Calibration **Analyst** Barnett, Brandi

**Due Date** // **Reason** **Complexity**

Scan the Laboratory Calibration Form and its attachments into the laboratory calibration folder in the image module. Using the edit findings function add a result to the evidence listed. In the results field enter the information for the simulator and simulator solution used to calibrate the instrument. The information listed should include: the name of the manufacturer of the simulator solution and level used, the lot number of the simulator solution, the simulator serial number and expiration date of the simulator calibration in the same format as below. Arial font, size 18, bold and centered.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The screenshot shows the 'General Analytical Module' window. The main area displays a tree view with the following structure:

- B-TEST
  - 0015 - Laboratory Calibration
    - 100666 - Test DMT - (B-TEST)
      - Answers-(Qty: 1)

Below the tree view, there are several controls:

- Result Type:** A dropdown menu set to 'Answers'.
- Instrument:** A dropdown menu with a small icon to its right.
- Qty:** A numeric input field set to '1'.
- Worklist ID:** A text input field with a red 'X' icon to its right.
- Font Settings:** A dropdown menu set to 'Arial' and a numeric input field set to '9'.
- Formatting Icons:** A row of icons for bold (B), italic (I), underline (U), text color (ABC), background color, bulleted list, numbered list, and indented list.
- Preview Icons:** Two icons for text alignment (A and A<sub>1</sub>).

The main display area shows the following text:

**Repco 0.100 Lot 9101**  
**Simulator 1234 Expires 2/2/2012**

At the bottom of the window, there are three buttons: 'Apply', 'Close', and an ellipsis (...). A small camera icon is located in the bottom right corner.

Press applies. This report is complete at administrative review. This request requires an administrative and technical review. The final report will be a certificate of laboratory calibration and will resemble the document below.



# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active



**State of Alaska**

Department of Public Safety

### SCIENTIFIC CRIME DETECTION LABORATORY REPORT

Sean Parnell, Governor  
Joseph A. Masters, Commissioner

#### Certificate of Calibration

##### Instrument Serial Number

100666

##### Simulator Solution and Simulator Information

**Repco 0.100 Lot 9101**

**Simulator 1234 Expires 2/2/2012**

The above listed instrument has been calibrated as of the date  
of this report per 13 AAC 63.100

A handwritten signature in cursive script, appearing to read "C. O'Byrne".

Forensic Scientist III: Brandi Barnett

Certification – This is the request used to certify an instrument for use after any necessary repairs are made. The request will already be made with the Check-In procedure; you will only need to assign it to yourself. Relate the evidence to the request. Scan the Certification Form and its attachments into the certification folder in the image module.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

**State of Alaska Scientific Crime Detection Laboratory**  
**DataMaster Certification Form – DMCEF2011R0**  
Issuing Authority – Nita Boiz, Scientific Director  
Approval date: 02/17/11 Effective date: 02/17/11 Archive date: Active

**DataMaster DMT Certification Form** Instrument # 100XXV  
Analyst Colleen  
Date 2/14/11

- Diagnostic Test (attached)
- Voltage Checks
  - Adjust voltages, if needed ☒
    - Lamp (1.01 V dc – 2.37 V dc)
    - Cooler (1.18 V dc – 2.74 V dc)
    - Bias (48 V dc – 112 V dc)
    - Chopper Frequency (475 hz – 575 hz)
  - Set RFI ☒
  - Adjust barometer reading, if needed ☒
  - Adjust volume setting, if needed ☒
  - Print Technician Screen (attached)
- Print Calibration Factors (attached)
- Volume Test - Minimum 1500 mL ☒
- Dry Gas Linearity Test (attached)
  - Levels Used ☒

	Mnf	Lot	Exp
i. 0.020g/210L	<u>Ilmo</u>	<u>1</u>	<u>2/14/12</u>
ii. 0.040g/210L	<u>Ilmo</u>	<u>2</u>	<u>2/14/12</u>
iii. 0.080g/210L	<u>Ilmo</u>	<u>3</u>	<u>2/14/12</u>
iv. 0.100g/210L	<u>Ilmo</u>	<u>4</u>	<u>2/14/12</u>
v. 0.200g/210L	<u>Ilmo</u>	<u>5</u>	<u>2/14/12</u>
vi. 0.300g/210L	<u>Ilmo</u>	<u>6</u>	<u>2/14/12</u>
- Functional Tests (attached)
  - Incomplete ☒
  - Invalid ☒
  - Interference ☒
    - Acetone - Mnf Repro Lot 9101 Level 0.040 Exp 2/12
      - Simulator Calibration Date 2/12 SN 1234
    - Isopropanol - Mnf Repro Lot 9101 Level 0.040 Exp 2/12
      - Simulator Calibration Date 2/12 SN 5678
  - RFI ☒
- Calibrate Stylus ☒
- Verify card reader is functional ☒
- Upload ☒
- Reset Options ☒
- Non-Drinking Subject Test (attached)
- NPAS Repair Paperwork (attached)

Using the edit findings function; add a result to the evidence listed. Then using the edit findings function, enter a value into findings for the evidence item. (A space in the results field will suffice as a result). No analytical result is needed. After you have entered the findings, press ok, it will ask if you want to create your cascading request of “Preventative Maintenance”, press yes and close the edit findings window. Press applies. This report is complete at administrative review. This request requires an administrative and technical review. The final report will be a certificate of laboratory certification and will resemble the document below.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active



### State of Alaska

Department of Public Safety

### SCIENTIFIC CRIME DETECTION LABORATORY REPORT

Sean Parnell, Governor  
Joseph A. Masters, Commissioner

February 24, 2011  
Requested Service: Certification

Laboratory Number: B-TEST



Serial #  
100666 DataMaster DMT

The accuracy of the calibration of the above listed instrument has been verified and therefore is certified for use as of the date of this report per 13 AAC 63.100.

A handwritten signature in cursive script, appearing to read "C. O'Bryant".

Forensic Scientist III: Brandi Barnett

Preventative Maintenance – This request is created automatically upon completion of a certification. The due date of the preventative maintenance request is set to four calendar years from the date of creation. If an instrument is brought back to the laboratory for preventative maintenance it should be checked in using normal check in procedures. The preventative maintenance request can be cancelled once the instrument has returned. A new request will be created upon completion of the certification.

# Alaska Scientific Crime Detection Laboratory


## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### ENTERING REPRESENTATIVES

Access the staff of the breath alcohol department using the administration tab, and choosing system administration.




Next chose staff from the laboratories options. Ensure the laboratory selection is “Breath” not “LAB”.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active



Press the green "+" at the bottom left of the staff box to add a new representative.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

**System Administration** **Laboratory Staff Listing**

List Maintenance

**Laboratories**

Breath

Barcodes

Department

Case Numbering

Event Notification

QueryTool DB Setup

Case Security Levels

Name	Title	Status
Abbott, Jacob Z	Operator	Active
Barnett, Brandi	Forensic Scientist III	Active
Bolz, Nita	Scientific Director	Active
Dym, Orin	Forensic Laboratory Manager	Active
Foster, Charles	Forensic Scientist III	Active
Nicholas, Dan	Forensic Technician	Active
O'Bryant, Colleen S	Forensic Scientist III	Active

Close ☒ Active only

Add <CTRL+A>

Filterable Labels

Activities

Archive Cases

Edit Evidence COC

Delete Cases

Delete Barcodes

Outside Agencies

Agency Info

Enter in the new staff member information using the key below

Archived 10

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

The screenshot shows a 'New Staff Member' form with the following fields and callouts:

- 1. Last name
- 2. First name
- 3. Middle initial
- 4. Title
- 5. Employee #
- 6. Status (Active)
- 7. Badge
- 8. Effective date (mm/dd/yyyy)
- 9. "Effective date" (chosed from drop down menu)
- 10. Expiration date (mm/dd/yyyy)
- 11. "Expiration date" (chosed from drop down menu)

Other fields include: E-mail, SSN, Supervisor, Mideo User Access Level, Time Zone Offset, Address 1, Address 2, City, State, ZIP, and a table for Phone Number, Type, and Description. The form has OK and Cancel buttons at the bottom.

1. Last name
2. First name
3. Middle initial
4. Operator/Supervisor
5. Agency
6. Active (always active)
7. Breath Alcohol ID#
8. Effective date (mm/dd/yyyy)
9. "Effective date" (chosed from drop down menu)
10. Expiration date (mm/dd/yyyy)
11. "Expiration date" (chosed from drop down menu)

Press OK to save the information.

# Alaska Scientific Crime Detection Laboratory

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

### APPENDIX III

#### BREATH ALCOHOL PROGRAM FORMS

##### Internal Forms

DataMaster Acceptance Form

Version: DMAF2011R0      Approved: 03/01/2011      Effective: 03/01/2011

DataMaster Check-In Form

Version: DMCIF2011R0      Approved: 03/01/2011      Effective: 03/01/2011

DataMaster Calibration Form

Version: DMCAF2011R0      Approved: 03/01/2011      Effective: 03/01/2011

DataMaster Certification Form

Version: DMCEF2011R0      Approved: 03/01/2011      Effective: 03/01/2011

##### External Forms (posted on the crime lab website)

Change in Instrument Status Form

Version: CIISF2011R0      Approved: 03/01/2011      Effective: 03/01/2011

Request for Operator Certification Form

Version: CIOSF2011R0      Approved: 03/01/2011      Effective: 03/01/2011



# Alaska Scientific Crime Detection Laboratory

---

## Breath Alcohol Procedure Manual

Issued: 3/01/2011  
Effective: 3/02/2011

Version: BAPM2011 R0  
Status: Active

Archived 10/17/2011