

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

MAY 12 2009

Supervisor/Operator Performing the Verification Procedure: DataMaster cdm S/N <u>130173</u>	
Name <u>J. Heinberg</u>	ID# <u>5828</u> Date <u>5-10-09</u>
A Agency <u>AST</u>	Phone # <u>822-3263</u>
Instrument Location <u>Glennallen</u>	
B Alco S/N <u>X172969</u> Target Value <u>0.078</u> High Pressure <u>1200</u>	
Alco Test Values <u>0.076</u> <u>0.076</u>	
1 st Alco 2 nd Alco	
Signature <u>[Signature]</u>	
(OVER) <u>[Signature]</u> <u>5/18/09</u>	

(Do Not write in the area below)


I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.

[Signature]
Nita J. Bolz
Scientific Director
State Breath Alcohol Program

Subscribed and sworn before me this 22 day of May, 2009.

[Signature] (Notary Seal Stamp)
Carolyn M. Noland
Notary Public, State of Alaska
Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N 130173 ✓

(CONTINUED FROM FRONT PAGE)

C TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130173

NOV 10, 2009 ✓ *COB*
SH/10

OPERATOR'S NAME:
SUPERVISOR *Heinbaugh*

OPERATOR'S NUMBER: 5828

SUBJECT'S LAST NAME:

UOC

SUBJECT'S FIRST NAME/MI :

UOC

O.L. #: UOC

DEPT/AGENCY: GAKO

CASE REPORT: UOC

TEST TYPE: U

ALCO TARGET VALUE: 0.078 ✓

ALCO 0-100 W172969 ✓

--- BREATH ANALYSIS ---

.078 ADJUSTED FOR 28.68 in ✓
ALCO TARGET .074 17:21
BLANK TEST .000 17:22
INTERNAL STANDARD VERIFIED 17:22
ALCO TU 28.68 in .076 ✓ 17:23
BLANK TEST .000 17:24
SUBJECT SAMPLE *COB* .000 17:24
SH/10 ✓
BLANK TEST .000 17:25
ALCO TU 28.68 in .076 17:25
BLANK TEST .000 17:26

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130173

NOV 10, 2009 ✓
TIME 17:29

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY

PROGRAM: OKAY

SOFTWARE DATE: 02/28/01

HEATERS

SAMPLE CHAMBER: 49c

BREATH TUBE: 41c

BAROMETER: 28.68 in

FLOW DETECTOR: OKAY

HIGH SPEED: OKAY

DETECTOR: OKAY

FILTERS: OKAY

QUARTZ STANDARD: OKAY

CALIBRATION: OKAY

PRINTER TEST

!"#\$%&'()*+,-./0123456789:;<=?@ABCDEFGHIJKL
MNPQRSTUVWXYZ[\]^_`abcdefgijklmnop
qrstuvwxyz{|}~■