

VERIFICATION OF CALIBRATION REPORT

MAY 28 2010

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Supervisor/Operator Performing the Verification Procedure:	
DataMaster cdm S/N <u>130177</u> ✓	
Name <u>STEVEN C. KEVAN</u>	ID# <u>6585</u> Date <u>5/24/10</u>
A Agency <u>ALASKA STATE TROOPERS</u> Phone # <u>438-2019</u>	
Instrument Location <u>ST. MARY'S</u> ✓	
B Alco S/N <u>X 301550</u> ✓ Target Value <u>.1079</u> ✓ High Pressure <u>500</u>	
Alco Test Value Average <u>.1081</u> ✓ <u>.1083</u> ✓	
1 st Alco 2 nd Alco	
Signature <u>[Signature]</u> ✓	
CAB 6/10/10	
(OVER)	

(Do Not write in the area below)

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.

[Signature]
Nita J. Bolz

Scientific Director
State Breath Alcohol Program

Subscribed and sworn before me this 24th day of June, 2010.

[Signature]
(Notary Seal Stamp)

Carolyn M. Noland
Notary Public, State of Alaska
Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N 130177 ✓

(CONTINUED FROM FRONT PAGE)

C TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

DIAGNOSTIC CHECK

Breath Test Record

Serial Number 130177 Model DataMaster
 Location SMD AST/St. Mary's CUB 6/10/10

Subject Information

Last Name Tst O.L. # 00
 First Name Tst

Operator Information

Operator Name Kevin/Steve Case 00
 Operator # 6585 Department SMYD

Test Results

Alco Target	<u>.079</u>	Alco Serial	<u>X301550</u>
Nominal	<u>0.079</u>	Barometer	<u>29.94</u>
Dry Gas 1	<u>0.081</u>	Time	<u>15:43:08</u>
Subject Alcohol	<u>0.000</u>	Time	<u>15:44:18</u>
Dry Gas 2	<u>0.083</u>	Time	<u>15:45:42</u>
		Report	<u>Y</u>

Diagnostic Test Record

Serial Number 130177 Model DataMaster
 Location SMD AST/St. Mary's

Diagnostic Test Results

Test Date	<u>05/24/2010</u>	Computer	<u>Okay</u>
Test Time	<u>15:38:27</u>	Program	<u>Okay</u>
		Software Date	<u>02/20/2001</u>
		Sample Chamber (C)	<u>49</u>
		Breath Tube (C)	<u>41</u>
		Barometer (in)	<u>29.94</u>
		Flow Detector	<u>Okay</u>
		High Speed Pump	<u>Okay</u>
		Detector	<u>Okay</u>
		Filter	<u>Okay</u>
		Quartz	<u>Okay</u>
		Calibration	<u>Okay</u>

- I, Nita J. Bolz
- (1) I am a F...
 - (2) The Ala...
 - (3) I am the...
 - (4) In that c...
 - responsible...
 - maintaining...
 - (5) The attac...
 - regularly rec...
 - verification.
 - (6) The refer...

VS:
 re De
 thin t
 ol pro
 raine
 ibrati
 oy a p
 tate c

.079 ADJUSTED FOR 29.94 in ✓
 ALCO TARGET .079 15:41
 BLANK TEST .000 15:42
 INTERNAL STANDARD VERIFIED 15:42
 ALCO TV 29.94 in .001 ✓ 15:43
 BLANK TEST .000 15:43
 SUBJECT SAMPLE .000 15:44
 BLANK TEST .000 15:45
 BLANK TEST .000 15:45
 BLANK TEST .000 15:46

Nita J. Bolz
 Scientific Director
 State Breath Alcohol Program

Subscribed and sworn before me this 08 day of June, 2010.

(Notary Seal Stamp)

Carolyn M. Noland
 Notary Public, State of Alaska
 Commission Expires with Office

FILTERS: OKAY
 QUARTZ STANDARD: OKAY
 CALIBRATION: OKAY
 PRINTER TEST
 !"#%&'()*+,-./0123456789:;<=>?@ABCDEF6
 HIJKLMNOPQRSTUVWXYZ[\]^_`abcdefg h i j k l m n o
 p q r s t u v w x y z { } ~