


# VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program


Supervisor/Operator Performing the Verification Procedure:		DataMaster cdm S/N	130178✓
Name	STEVE DUNN	ID#	4106✓
Date	12-8-10✓		
A Agency	APD	Phone #	786-8900
Instrument Location	EAGLE RIVER SUB		
B Alco S/N	75260✓	Target Value	.081✓
High Pressure	700 PS		
Alco Test Values	.082✓	.082✓	
1st Alco		2nd Alco	
Signature			
(OVER)			

(Do Not write in the area below)

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:


- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.

  
Nita J. Bolz  
Scientific Director  
State Breath Alcohol Program

Subscribed and sworn before me this 18th day of JAN, 2010. 

Daniel L. Nicholas  
Notary Public, State of Alaska  
Commission Expires with Office

(Notary Seal Stamp)



# VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N

130178

8mb  
12/23/10

(CONTINUED FROM FRONT PAGE)

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

## NONDRINKING SUBJECT TEST

## DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY  
DATAMASTER cdm 130178

DECEMBER 08, 2010

OPERATOR'S NAME:

DUNN/S

OPERATOR'S NUMBER: 4106

SUBJECT'S LAST NAME:

VERIFICATION OF CAL

SUBJECT'S FIRST NAME/MI :

V

O.L. #: V

DEPT/AGENCY: ANCL

CASE/REPORT: V

TEST TYPE: V

ALCO TARGET VALUE: .081

ALCO S/N: 75260

--- BREATH ANALYSIS ---

.081 ADJUSTED FOR 29.47 in  
ALCO TARGET .079  
BLANK TEST .000  
INTERNAL STANDARD VERIFIED  
ALCO TV 29.47 in .082  
BLANK TEST .000  
SUBJECT SAMPLE .000  
BLANK TEST .000  
ALCO TV 29.47 in .082  
BLANK TEST .000

ALASKA DEPARTMENT OF PUBLIC SAFETY  
DATAMASTER cdm 130178

DECEMBER 08, 2010  
TIME 13:14

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY  
PROGRAM: OKAY  
SOFTWARE DATE: 01/08/09  
HEATERS  
SAMPLE CHAMBER: 49c  
BREATH TUBE: 42c  
BAROMETER: 29.49 in  
FLOW DETECTOR: OKAY  
PUMP  
HIGH SPEED: OKAY  
DETECTOR: OKAY  
FILTERS: OKAY  
QUARTZ STANDARD: OKAY  
CALIBRATION: OKAY

PRINTER TEST

!"#\$%&'()\*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNO  
PQRSTUVWXYZ[\]^\_`abcdefghijklmnopqrstuvwxyz{|}~Δ