

OCT 31 2011

VERIFICATION OF CALIBRATION REPORT OCT 31 2011

*of DataMaster cdm Breath Test Instrument
State of Alaska*

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N <u>130 281</u>										
Supervisor/Operator Performing the Verification Procedure:										
Name <u>Trapper James Kimura</u>	ID# <u>5677</u> Date <u>10/28/11</u>									
A Agency <u>Alaska State Troopers</u>	Phone # <u>225-5118</u>									
Instrument Location <u>Ketchikan AST Post</u>										
B Alco S/N <u>78268</u> Target Value <u>0.081</u> High Pressure <u>300 PSI</u>										
<table style="width: 100%; border: none;"> <tr> <td style="width: 30%;"></td> <td style="width: 30%; text-align: center;">0.079 <small>Corr 10/01/11</small></td> <td style="width: 30%;"></td> </tr> <tr> <td>Alco Test Values</td> <td style="text-align: center;"><u>.081</u></td> <td style="text-align: center;"><u>.081</u></td> </tr> <tr> <td></td> <td style="text-align: center;"><small>1st Alco</small></td> <td style="text-align: center;"><small>2nd Alco</small></td> </tr> </table>			0.079 <small>Corr 10/01/11</small>		Alco Test Values	<u>.081</u>	<u>.081</u>		<small>1st Alco</small>	<small>2nd Alco</small>
	0.079 <small>Corr 10/01/11</small>									
Alco Test Values	<u>.081</u>	<u>.081</u>								
	<small>1st Alco</small>	<small>2nd Alco</small>								
Signature <u>[Handwritten Signature]</u>	<u>[Handwritten Signature]</u> 10/28/11									

(OVER)

(Do Not write in the area below)


I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.

[Handwritten Signature]
 Nita J. Bolz
 Scientific Director
 State Breath Alcohol Program

Subscribed and sworn before me this 7 day of Nov, 2011 AK

[Handwritten Signature] (Notary Seal Stamp)
 Daniel L. Nicholas Nikki Roth
 Notary Public, State of Alaska
 Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument
State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N 130289

COB 10/31/11

(CONTINUED FROM FRONT PAGE)

C TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130289

OCTOBER 28, 2011 COB 10/31/11

OPERATOR'S NAME:
KINURA/JAMES/K
OPERATOR'S NUMBER: 5677
SUBJECT'S LAST NAME:
VOC
SUBJECT'S FIRST NAME/MI :
VOC
O.L. #: 000000
DEPT/AGENCY: KETA
CASE/REPORT: 11-000000
TEST TYPE: V
ALCO TARGET VALUE: .061
ALCO S/N: 76268

--- BREATH ANALYSIS ---

.061 ADJUSTED FOR 29.38 in		
ALCO TARGET	.079	10:01
BLANK TEST	.000	10:02
INTERNAL STANDARD	VERIFIED	10:02
ALCO TV 29.38 in	.061	10:03
BLANK TEST	.000	10:03
SUBJECT SAMPLE COB 10/31/11	.060	10:04
BLANK TEST	.000	10:04
ALCO TV 29.39 in	.061	10:05
BLANK TEST	.000	10:05

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130289

OCTOBER 28, 2011
TIME 10:07

--- DIAGNOSTIC CHECK ---

COMPUTER:	OKAY
PROGRAM:	OKAY
SOFTWARE DATE:	01/08/09
HEATERS	
SAMPLE CHAMBER:	45c
BREATH TUBE:	41c
BAROMETER:	29.39 in
FLOW DETECTOR:	OKAY
PUMP	
HIGH SPEED:	OKAY
DETECTOR:	OKAY
FILTERS:	OKAY
QUARTZ STANDARD:	OKAY
CALIBRATION:	OKAY

PRINTER TEST

!@#x"()+-./0123456789:;<=>?@BCDEFG
HIJKLMNPQRSTVWXYZ[-_`abcde fghijklmno
pqrsuvwxyz{|}~"

