

OCT 25 2010

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska


Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Supervisor/Operator Performing the Verification Procedure:		DataMaster cdm S/N	130292 ✓
Name	TRACY MORLEY ✓	ID#	7252 ✓
Date	10-15-10 ✓		
A	Agency	JUNEAU Police Dept	Phone #
	Instrument Location	JPD - JUNEAU	
B	Alco S/N	X 172933 ✓	Target Value
		.078 ✓	High Pressure
		30.20	
	Alco Test Values	.079 ✓	.080 ✓
		1st Alco	2nd Alco
	Signature	COB 11/9/10	
		(OVER)	


(Do Not write in the area below)


I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.


Nita J. Bolz
Scientific Director
State Breath Alcohol Program

Subscribed and sworn before me this 9th day of NOV, 2010.


(Notary Seal Stamp)
Daniel L. Nicholas
Notary Public, State of Alaska
Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N

130292

(CONTINUED FROM FRONT PAGE)

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130292

OCTOBER 15, 2010 ^{DOB} 11/9/10

OPERATOR'S NAME:
MURPHY/TRACY/ERNEST
OPERATOR'S NUMBER: 7252
SUBJECT'S LAST NAME:
VO
SUBJECT'S FIRST NAME/MI :
VOC
O.L. #: VOC
DEPT/AGENCY: JNU1
CASE/REPORT: VOC
TEST TYPE: V
ALCO TARGET VALUE: .078
ALCO S/N: X172933

--- BREATH ANALYSIS ---

.078 ADJUSTED FOR 30.20 in		
ALCO TARGET	.078	21:10
BLANK TEST	.000	21:12
INTERNAL STANDARD	VERIFIED	21:12
ALCO TV 30.20 in	.079	21:12
BLANK TEST	.000	21:13
SUBJECT SAMPLE	.000	21:13
BLANK TEST	.000	21:14
ALCO TV 30.20 in	.080	21:14
BLANK TEST	.000	21:15

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130292

OCTOBER 15, 2010
TIME 21:16

--- DIAGNOSTIC CHECK ---

COMPUTER:	OKAY
PROGRAM:	OKAY
SOFTWARE DATE:	01/08/09
HEATERS	
SAMPLE CHAMBER:	48c
BREATH TUBE:	41c
BAROMETER:	30.20 in
FLOW DETECTOR:	OKAY
PUMP	
HIGH SPEED:	OKAY
DETECTOR:	OKAY
FILTERS:	OKAY
QUARTZ STANDARD:	OKAY
CALIBRATION:	OKAY

PRINTER TEST

10122455789...@ARCT

noticed
printer test
11/9/10