

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N

~~130256~~ ^{02A} 9/3/10
130526

(CONTINUED FROM FRONT PAGE)

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130525

AUGUST 30, 2010 ^{COB 9/7/10}

OPERATOR'S NAME:
SUPERVISOR
OPERATOR'S NUMBER: 4188
SUBJECT'S LAST NAME:
VERIFICATION OF CAL
SUBJECT'S FIRST NAME/MI :
L
O.L. #: L
DEPT/AGENCY: SDP1
CASE/REPORT: K
TEST TYPE: V
ALCO TARGET VALUE: .079
ALCO S/N: K173004

--- BREATH ANALYSIS ---

| | | |
|----------------------------|----------|-------|
| .079 ADJUSTED FOR 30.20 in | | |
| ALCO TARGET | .079 | 17:14 |
| BLANK TEST | .000 | 17:16 |
| INTERNAL STANDARD | VERIFIED | 17:16 |
| ALCO TV 30.20 in | .077 | 17:16 |
| BLANK TEST | .000 | 17:17 |
| SUBJECT SAMPLE | .000 | 17:17 |
| BLANK TEST | .000 | 17:18 |
| ALCO TV 30.20 in | .078 | 17:19 |
| BLANK TEST | .000 | 17:19 |

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130525

AUGUST 30, 2010
TIME 17:21

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY
PROGRAM: OKAY
SOFTWARE DATE: 02/20/01
HEATERS
SAMPLE CHAMBER: 50c
BREATH TUBE: 41c
BAROMETER: 30.20 in
FLOW DETECTOR: OKAY
PUMP
HIGH SPEED: OKAY
DETECTOR: OKAY
FILTERS: OKAY
QUARTZ STANDARD: OKAY
CALIBRATION: OKAY

PRINTER TEST

! " # \$ % & ' () * + , - . / 0 1 2 3 4 5 6 7 8 9 : ; < = > ? @ A B C D E F G
H I J K L M N O P Q R S T U V W X Y Z I \] ^ _ ` a b c d e f g h i j k l m n o
p q r s t u v w x y z { | } ~