

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

MAR 23 2010

Supervisor/Operator Performing the Verification Procedure: DataMaster cdm S/N 130666 ✓

Name John Merrion ID# 4785 Date 3-17-10

A Agency Fairbanks Police Department Phone # 450-6500

Instrument Location 911 Cushman Street Fairbanks, Alaska 99701

B Alco S/N X172971 ✓ Target Value .083 ✓ High Pressure 800

Alco Test Value Average .084 ✓ .085 ✓
1st Alco 2nd Alco

Signature [Signature] ✓

(OVER)

AMB 4/1/10

(Do Not write in the area below)

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.

[Signature]
Nita J. Bolz
Scientific Director
State Breath Alcohol Program

Subscribed and sworn before me this 5th day of April, 2010.

[Signature]
Carolyn M. Noland
Notary Public, State of Alaska
Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N 130666 ✓

(CONTINUED FROM FRONT PAGE)

C TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130666 ✓

MARCH 17, 2010

OPERATOR'S NAME:
MERRION/JOHN/E
OPERATOR'S NUMBER: 4785
SUBJECT'S LAST NAME:
VERIFICATION/FOCAL
SUBJECT'S FIRST NAME/MI:
A

O.L. #: A
DEPT/AGENCY: FBK1
CASE/REPORT:
TEST TYPE: V
ALCO TARGET VALUE: .083 ✓
ALCO S/N: X172971 ✓

--- BREATH ANALYSIS ---

.083 ADJUSTED FOR 30.25 in ✓
ALCO TARGET .083 21:23
BLANK TEST .000 21:23
INTERNAL STANDARD VERIFIED 21:24
ALCO TV 30.25 in .084 21:24
BLANK TEST .000 21:25
SUBJECT SAMPLE .000 21:25
BLANK TEST .000 21:26
ALCO TV 30.25 in .085 21:26
BLANK TEST .000 21:27

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130666 ✓

MARCH 17, 2010 ✓
TIME 21:19

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY
PROGRAM: OKAY
SOFTWARE DATE: 02/20/01
HEATERS
SAMPLE CHAMBER: 43c
BREATH TUBE: 43c
BAROMETER: 30.25 in
FLOW DETECTOR: OKAY
PUMP
HIGH SPEED: OKAY
DETECTOR: OKAY
FILTERS: OKAY
QUARTZ STANDARD: OKAY
CALIBRATION: OKAY

PRINTER TEST

!@#%&'()*+,-./0123456789:;<=>?@BCDEFG
HIJKL MNOPQRSTU VWXYZ[\]^_`abcde fghijklmno
pqrstuvwxyz{|}~"