## VERIFICATION OF CALIBRATION REPORT

SEP 1 4 2009

of DataMaster cdm Breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Supervisor/Operator Performing the Verification Procedure:  Name Doreen Valadez  Agency 6eWard PD  DataMaster cdm S/N /30/3  ID# 3423 Date 9/	409
Instrument Location 410 Adams () 1111	3338
Alco S/N X 301605 Target Value 1082 High Pressure 300	
Alco Test Value Average  Alco Que Alco	<u>psi</u>
Signature	
(OVER)	10/3/09

(Do Not write in the area below)

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.

(5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Scientific Director


State Breath Alcohol Program

Subscribed and sworn before me the

Notary Seal Stamp)

Carolvn M. Noland

Notary Public, State of Alaska Commission Expires with Office



## VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N /3()/32

(CONTINUED FROM FRONT PAGE)
Datamaster moved upstairs for two weeks during jail repairs.

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

## NONDRINKING SUBJECT TEST DIAGNOSTIC CHECK ALASKA DEPARTMENT OF PUBLIC SAFETY DATAMASTER cdm 130132 ALASKA DEPARTMENT OF PUBLIC SAFETY SEPTEMBER 10, 2009 DATAMASTER cdm 130132 TIME 12:47 SEPTEMBER 10, 2009 --- DIAGNOSTIC CHECK ---OPERATOR'S NAME: COMPUTER: OKRY XXX OPERATOR'S NUMBER: 3423 SUBJECT'S LAST NAME: VERIFICATION OF CAL PROGRAM: OKRY SOFTWARE DATE: 03/29/00 SUBJECT'S FIRST NAME/MI: HEATERS O.L. #: X DEPT/AGENCY: SMD1 Sample Chamber: BREATH TUBE: CASE/REPORT: 0 TEST TYPE: V BAROMETER: 29.15 in ALCO TARGET VALUE: .082 HTCO 274: X201602 FLOW DETECTOR: THE STATE OF BREATH ANALYSIS PIMP HIGH SPEED: OKAY .082 ADJUSTED FOR 29.15 in 12:39 ALCO TARGET DETECTOR: OKRY BLANK TEST .966 12:40 12:49 INTERNAL STANDARD VERIFIED) FILTERS: OKRY ALCO TV 29.15 in .977 12:40 BLANK TEST .000 12:41 QUARTZ STANDARD: OKRY SUBJECT SAMPLE .000 12:43 BLANK TEST .000 12:44 CALIBRATION: OKAY ALCO TV 29.15 in 12:44 .876 12:45 **BLANK TEST** .000 PRINTER TEST !"#\$%%'()\*+,-./0123456789:;<=>?@ABCDEFG HIJKLMNOPORSTUVNKYZ[\]^\_\*abcdefghijklmno parsturexyz(!)~