VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument State of Alaska

FEB 0 8 2010

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Scientific Crin	de Detection Laboratory - Statewide Breath Aconst Program
Supervisor/Operator Performing	
Name STEVE	Phone # 800-478-4112
Agency	BETHEL AST
Alco S/N 7767	Target Value
Alco Test Value Averag	$\frac{.09\lambda}{1^{st} \text{ Alco}} \frac{.091}{2^{nd} \text{ Alco}}$
Signature	Stock
	(OVER)
	(Do Not write in the area below)
 (2) The Alaska Scientific Crime Detection (3) I am the Scientific Director of the State (4) In that capacity, I am responsible for or responsible for verifying the calibration of maintaining the records of the program. (5) The attached verification is a true and regularly recorded activity of the breath a 	e of Alaska Scientific Crime Detection Laboratory. Laboratory is an entity within the Department of Public Safety.
verification. (6) The referenced instrument is certified Nita J. Bolz	for evidentiary use in the State of Alaska.
Scientific Director State Breath Alcohol Program Subscribed and sworn before me this	day of 60., 2010.
1 \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \	Seal Stamp)
Notary Public, State of Alaska Commission Expires with Office	

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N _ 130 13Le

(CONTINUED FROM FRONT PAGE)

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

DIAGNOSTIC CHECK



ALASKA DEPARTMENT (DATAMASTER cde 130)	
FEBRUARY 05, 2010 TIME 11:58	
DIAGNOST	TIC CHECK
COMPUTER:	OKRY
PROGRAM:	OKAY
SUFTWARE DATE:	01/08/09
HEATERS SAMPLE CHAMBER: BREATH TUBE:	49c 48c
BAROMETER:	29.73 in
FLOW DETECTOR:	OKAY
PUMP HIGH SPEED:	OKAY
DETECTOR:	OKAY
FILTERS:	OKAY
QUARTZ STANDARD:	OKAY
CALIBRATION:	OKAY
PRINTE !"#\$%%?()++,-,20123 HIJKLMNOPGRSTUVUXYZC pgrstuvuxyz({}}~#	R TEST 456789::<=>?@ABCDEF \J^_'abcdefghijklmn