## VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument State of Alaska

SEP 10 200

	visor/Operator Performing the Verification Procedure:  Name
A	Agency Phone # <u>786 - 8 9</u>
	Instrument Location
В	Alco S/N X 30/563 Target Value . 079 High Pressure 500 P
	Alco Test Value Average OSO 2nd Alco
	Signature
	(OVER)
lz afte	(Do Not write in the area below) or being first duly sworn, depose and state as follows:

(5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the

verification.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Wita J. Bolz Scientific Director

State Breath Alcohol Program

Subscribed and sworn before me this

(Notary Seal Stamp)

Serrolyn M. Noland

Notary Public, State of Alaska Commission Expires with Office

## VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

(CONTINUED FROM FRONT PAGE)

 $\mathbf{C}$ 

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

## NOND, INVING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY DATAMASTER cdm 130579

SEPTEMBER 09, 2009

OPERATOR'S NAME:

DUNN/S

OPERATOR'S NUMBER: 4106

SUBJECT'S LAST NAME:

VERIFICATION OF CAL

SUBJECT'S FIRST NAME/MI :

V

O.L. #: V

DEPT/AGENCY: ANC1

CASE/REPORT: V

TEST TYPE: V

ALCO TARGET VALUE: .079 ~

ALCO S/N: X301563 🗸

--- BREATH ANALYSIS ---

.079 ADJUSTED FOR 29.67 in

ALCO TARGET .078 23:27 BLANK TEST .000 23:28 INTERNAL STANDARD 23:28

VERIFIED ALCO TV 29.67 in .080 🗸

23:29 BLANK TEST .000 23:29

SUBJECT SAMPLE 4 .000 23:30 BLANK TEST .000 23:31

ALCO TV 29.68 in .079 V 23:31 BLANK TEST .000 23:32 DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY DATAMASTER cdm 130579

SEPTEMBER 09, 2009 🗸 TIME 23:23

--- DIAGNOSTIC CHECK ---

COMPUTER:

OKAY

PROGRAM:

OKAY

SOFTWARE DATE:

12/11/08

**HEATERS** 

SAMPLE CHAMBER:

50c

BREATH TUBE:

42c

BAROMETER:

29.67 in

FLOW DETECTOR:

OKAY

PUMP

HIGH SPEED:

OKAY

DETECTOR:

OKAY

FILTERS:

OKAY

OKAY

CALIBRATION:

QUARTZ STANDARD:

OKAY

PRINTER TEST

!"#\$%&'()\*+,-./0123456789:;<=>?@ABCDEFG HIJKLMNOPQRSTUVWXYZ[\]^\_ abcdefghijklmno pgrstuvwxyz{|}~\dagger

S D-+1361