

# VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

AUG 31 2009

|                                                                                    |                                      |
|------------------------------------------------------------------------------------|--------------------------------------|
| Supervisor/Operator Performing the Verification Procedure:                         |                                      |
| DataMaster cdm S/N <u>130196</u> ✓                                                 |                                      |
| Name <u>J Swartz</u>                                                               | ID# <u>4314</u> Date <u>03/30/09</u> |
| A Agency <u>SCDL</u>                                                               | Phone # <u>269-5592</u>              |
| Instrument Location <u>SCDL</u> ✓                                                  |                                      |
| B Alco S/N <u>SC0108</u> ✓ Target Value <u>.080</u> ✓ High Pressure <u>850 PSI</u> |                                      |
| Alco Test Values <u>.080</u> ✓ / <u>.080</u> ✓                                     |                                      |
| 1 <sup>st</sup> Alco / 2 <sup>nd</sup> Alco                                        |                                      |
| Signature <u>J Swartz</u> ✓                                                        |                                      |
| (OVER)                                                                             |                                      |

BMB  
9/10/09


I, Nita J. Bolz, after being first duly sworn, depose and state as follows (Do not write in the area below)

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.
- (6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Nita J. Bolz  
Nita J. Bolz  
Scientific Director  
State Breath Alcohol Program

Subscribed and sworn before me this 29<sup>th</sup> day of Sept, 2009.

Carolyn M. Noland (Notary Seal Stamp)  
Carolyn M. Noland  
Notary Public, State of Alaska  
Commission Expires with Office



# VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N 130194 —

(CONTINUED FROM FRONT PAGE)

C TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

## NONDRINKING SUBJECT TEST

ALASKA DEPARTMENT OF PUBLIC SAFETY  
DATAMASTER cdm 130194

MARCH 30, 2009

OPERATOR'S NAME:  
SWARTZ/J

OPERATOR'S NUMBER: 4314

SUBJECT'S LAST NAME:

VERIFICATION OF CAL

SUBJECT'S FIRST NAME/NI :

L

O.L. #: L

DEPT/AGENCY: AA01

CASE/REPORT:

TEST TYPE: V

ALCO TARGET VALUE: .000

ALCO S/N: SCOTT08

### — BREATH ANALYSIS —

.000 ADJUSTED FOR 29.44 in

ALCO TARGET .078 09:27

BLANK TEST .000 09:27

INTERNAL STANDARD VERIFIED 09:28

ALCO TV 29.44 in .000 09:28

BLANK TEST .000 09:29

SUBJECT SAMPLE .000 09:30

BLANK TEST .000 09:31

ALCO TV 29.44 in .000 09:31

BLANK TEST .000 09:32

8/3/20/09

## DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY  
DATAMASTER cdm 130194

MARCH 30, 2009

TIME 09:25

### — DIAGNOSTIC CHECK —

COMPUTER: OKAY

PROGRAM: OKAY

SOFTWARE DATE: 01/08/09

HEATERS

SAMPLE CHAMBER: 40c

BREATH TUBE: 41c

BAROMETER: 29.44 in

FLOW DETECTOR: OKAY

PUMP  
HIGH SPEED: OKAY

DETECTOR: OKAY

FILTERS: OKAY

QUARTZ STANDARD: OKAY

CALIBRATION: OKAY

### PRINTER TEST

!"#\$%&'()\*+,-./0123456789:;<=?@ABCDEFGHIJ  
KLMNOPQRSTUVWXYZ[\]^\_`abcdefghijklmnopqr  
stuvwxyz{|}~