

Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form

Effective: 1/18/2024

Version: 4.0

Instrument #: 100418

Analyst/Technician(s): Carlie Blaister

Start Date: 12/23/25

Date Instrument Received: 12/23/25

Reason for instrument evaluation

Screen issues (frozen), not corrected w/ reboots.

Visual Inspection:

☒ Clean

☒ Replace missing screws

☒ Cermetex Box (modem) is secure

☒ Ensure modem cable present

☒ Double Check Valve secured to wall

☒ All tubing and connections secure

☒ O-ring inspection and ensure 2 replacement O-rings in bag

☐ Attach VOC

☐ VOC Internal Standard $\leq 3.00\%$

☐ Verify software version

Condition of item: ☐ No apparent damage ☒ Damaged

OK 12/23/25

NOTE: As Found Linearity is required if the instrument is not damaged and the internal standard value is not more than 4.00%

As Found Conditions (if performed):

Barometric Pressure (inHg): _____ Temperature (C): _____ Relative Humidity (%) _____

OK 12/23/25

Linearity Test Levels (g/210L): 0.020, 0.080, 0.300

Notes: Keyboard works to wake screen. Touch screen not working. Unable to use instrument

☐ Attach Non-Drinking Subject Test

OK 12/23/25

Conclusion:

☐ Ready for Adjustment

☒ Needs Repair

☐ Ready for Calibration

Lab Work End Date: 12/23/25

Page 1 of 1

All printed copies are uncontrolled.

Approved by: Chemistry Supervisor