

# Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form

Effective: 3/20/2025

Version: 4.0

Instrument #: 100381

Received From: APD

Date: 6/23/25

Analyst: C. O. Bryant

Reason for instrument return:

Drifting detector voltage on graph during low  
or no breath flow.

### From the field

Change instrument and case location in JusticeTrax Breath Alcohol Database

Ensure change in instrument status form is in JusticeTrax

Upload

Update agency information and polling status in DMHost Calling Program

Check for VOCs

### From Intoximeters

Intoximeters Repair Paperwork (attached)

Change instrument and case location in JusticeTrax Breath Alcohol Database

Reset Options and Clear Memory