

# Alaska Scientific Crime Detection Laboratory

## Change in Instrument Status Form

Issued: 12/15/2017  
Effective: 12/15/2017

Version: CIISF 2017 R0  
Status: Active

Date: 02/11/2020 Instrument Serial Number: 100378

Instrument Location: AST/Fairbanks

Supervisor Name and Agency: Colleen O'Bryant, SCDL

Describe the Status Change or Issue with the Instrument:

---

---

---

---

---

---

---

---

---

---

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: NA

Out of Service Date/Time: NA

-----  
**For Use by SCDL**

**Additional Notes**

See the attached email documenting gas flow and pump errors during the 1/31/2020 - 02/02/2020 time frame as well as the troubleshooting performed by Trooper David Bozman (Breath Test Supervisor). I confirmed that following the troubleshooting there have not been any additional occurrences of gas flow or pump errors.

COB 02/11/2020

**Email completed form to [dps.scdl.toxicology@alaska.gov](mailto:dps.scdl.toxicology@alaska.gov)**

For questions contact the Breath Alcohol Program at 907-269-5740

## O'Bryant, Colleen S (DPS)

---

**From:** Barnett, Brandi M (DPS)  
**Sent:** Monday, February 03, 2020 9:48 AM  
**To:** O'Bryant, Colleen S (DPS)  
**Subject:** FW: Fairbanks AST Datamaster

Can you follow up on this and make sure we get some form of documentation?

Thanks,  
Brandi

**From:** Bozman, David W (DPS) <david.bozman@alaska.gov>  
**Sent:** Monday, February 3, 2020 9:39 AM  
**To:** Barnett, Brandi M (DPS) <brandi.barnett@alaska.gov>  
**Subject:** Fairbanks AST Datamaster

Good Morning Brandi,

Over the last week we have gotten a Flow Error and a Pump Error on the Fairbanks AST Datamaster. I looked at it this morning and the check valve had fallen off the side of the unit. I fixed it with double sided tape and ran a test run after fixing it. I will let you know if there are any other issues.

Trooper David Bozman  
Fairbanks Patrol  
"D" Detachment  
1979 Peger Road  
Fairbanks, AK99709  
Office: (907)-451-5100  
Cell: (907)-953-8015  
Fax: (907)451-5165

Serial Number	Test Date	Location	Status
100378	1/1/2020 9:05 AM	AST Fairbanks	
100378	1/1/2020 9:36 PM	AST Fairbanks	
100378	1/1/2020 10:16 ...	AST Fairbanks	
100378	1/4/2020 1:46 AM	AST Fairbanks	
100378	1/4/2020 1:53 AM	AST Fairbanks	
100378	1/7/2020 10:55 ...	AST Fairbanks	
100378	1/15/2020 8:26 ...	AST Fairbanks	RFI
100378	1/15/2020 8:32 ...	AST Fairbanks	RFI
100378	1/15/2020 8:35 ...	AST Fairbanks	
100378	1/16/2020 1:59 ...	AST Fairbanks	
100378	1/16/2020 2:11 ...	AST Fairbanks	Maximum Attempts Exceeded
100378	1/16/2020 2:17 ...	AST Fairbanks	RFI
100378	1/17/2020 7:32 ...	AST Fairbanks	
100378	1/19/2020 12:24...	AST Fairbanks	
100378	1/19/2020 4:47 ...	AST Fairbanks	
100378	1/19/2020 5:04 ...	AST Fairbanks	Maximum Attempts Exceeded
100378	1/19/2020 5:11 ...	AST Fairbanks	Maximum Attempts Exceeded
100378	1/21/2020 9:52 ...	AST Fairbanks	
100378	1/21/2020 10:02...	AST Fairbanks	
100378	1/25/2020 9:52 ...	AST Fairbanks	
100378	1/26/2020 2:49 ...	AST Fairbanks	
100378	1/26/2020 4:50 ...	AST Fairbanks	RFI
100378	1/26/2020 4:54 ...	AST Fairbanks	
100378	1/31/2020 11:28...	AST Fairbanks	Gas Flow Error
100378	1/31/2020 11:30...	AST Fairbanks	Gas Flow Error
100378	1/31/2020 11:38...	AST Fairbanks	Gas Flow Error
100378	2/2/2020 9:56 PM	AST Fairbanks	Pump Error
100378	2/2/2020 9:58 PM	AST Fairbanks	Gas Flow Error
100378	2/2/2020 10:01 ...	AST Fairbanks	Gas Flow Error
100378	2/2/2020 10:05 ...	AST Fairbanks	Gas Flow Error
100378	2/3/2020 9:14 AM	AST Fairbanks	
100378	2/3/2020 9:27 AM	AST Fairbanks	
100378	2/3/2020 11:31 ...	AST Fairbanks	
100378	2/3/2020 9:12 PM	AST Fairbanks	
100378	2/3/2020 9:23 PM	AST Fairbanks	
100378	2/6/2020 5:29 PM	AST Fairbanks	
100378	2/8/2020 4:42 AM	AST Fairbanks	
100378	2/8/2020 5:39 AM	AST Fairbanks	RFI
100378	2/8/2020 5:40 AM	AST Fairbanks	
100378	2/9/2020 4:56 PM	AST Fairbanks	
100378	2/9/2020 5:06 PM	AST Fairbanks	