

Alaska Scientific Crime Detection Laboratory

Change in Instrument Status Form

Issued: 6/11/2014
Effective: 6/11/2014

Version: CIISF 2014 R0
Status: Active

Date: 11/28/16 Instrument Serial Number: 100386

Instrument Location: UAF PD

Supervisor Name and Agency: Sam Beaudreault, UAF PD
(see email)

Describe the Status Change or Issue with the Instrument:

Software/USB issue originating with the
card reader.

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: _____

Out of Service Date/Time: N/A

For Use by SCDL

Additional Notes

Restarted the instrument and completed a
test without issue. CBS 11/28/16

Email completed form to dps.scdl.toxicology@alaska.gov

For questions contact the Breath Alcohol Program at 907-269-5740

O'Bryant, Colleen S (DPS)

From: DPS SCDL Toxicology (DPS sponsored)
Sent: Monday, November 28, 2016 7:15 AM
To: Foster, Charles R (DPS); O'Bryant, Colleen S (DPS); Barnett, Brandi M (DPS)
Subject: FW: Datmaster glitch

From: Samuel Beaudreault [sabeaudreault@alaska.edu]
Sent: Thursday, November 24, 2016 2:00 PM
To: DPS SCDL Toxicology (DPS sponsored)
Subject: Datmaster glitch

On 11/23/16 at about 3 am our Datmaster dmt malfunctioned.

AST was borrowed our (UAFPD's) Datmaster.

Everything was going fine until he used the card scan port on the left side of the machine for the subjects driver's license.

That is when the screen froze and did not allow anymore input, no matter what he or our officer did to it.

They got it working again by turning the power off and on, and it worked fine.

I did a test run just now and it seems to be working correctly.

- FYI

and I hope you had a good Thanksgiving!

--
Officer Sam Beaudreault
UAF Police Department
ph. (907)474-7721
fax (907)474-5555