

Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form

Effective: 11/10/2021

Version: 2.0

Instrument #: 100683

Analyst(s): Carlisle Glaister

Start Date: 7/11/2022

Reason for instrument evaluation

Screen Issues

Visual Inspection

- | | |
|---|---|
| <input checked="" type="checkbox"/> Clean | <input checked="" type="checkbox"/> Replace missing screws |
| <input checked="" type="checkbox"/> Cermetex Box (modem) is secure | <input checked="" type="checkbox"/> Ensure modem cable present |
| <input checked="" type="checkbox"/> Double Check Valve secured to wall | <input checked="" type="checkbox"/> All tubing and connections secure |
| <input checked="" type="checkbox"/> O-ring inspection and ensure 2 replacement O-rings in bag | |

Troubleshooting/Repair

Screen not responding to touch. Recommend screen replacement

Unable to run VOC or Non-Drinking Subject Test

- VOC (attached) Verify software version ENG 7/11/2022
 VOC Internal Standard value $\leq 3.00\%$ Non-Drinking Subject Test (attached)

- | | |
|---|--|
| <input type="checkbox"/> Return to Factory | <input type="checkbox"/> Ready for Calibration |
| <input type="checkbox"/> Ready for Certification | <input type="checkbox"/> Proceed to Close Out |
| <input checked="" type="checkbox"/> Further Evaluation Needed | |

End Date: 7/11/2022

Page 1 of 1

All printed copies are uncontrolled.

Approved by: Chemistry Supervisor