

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Issued: 1/13/2017
Effective: 1/17/2017

Version: DMCIF 2017 R0
Status: Active

Instrument S/N 100373

Received From APD Training

Date 04/26/21

Analyst Derch Walton

1. Reason for instrument return:

Returned from training

2. Change instrument and case location in JusticeTrax Breath Alcohol Database

3. Ensure change in instrument status form is in JusticeTrax

4. Upload

Instrument was not in the field. Reset options and cleared memory. JHW 04/26/21

5. Update agency information and polling status in DM Host Calling Program

6. Check for VOCs

From Intoximeters

7. Intoximeters Repair Paperwork (attached)

Instrument needs to be returned to the manufacturer for refurbishment. JHW 04/26/21