

Alaska Scientific Crime Detection Laboratory

Change in Instrument Status Form

Issued: 12/15/2017
Effective: 12/15/2017

Version: CIISF 2017 R0
Status: Active

Date: 6/29/20 Instrument Serial Number: 100406

Instrument Location: Juneau Police Dept

Supervisor Name and Agency: Sterling Salisbury, JPD

Describe the Status Change or Issue with the Instrument:

I was advised by nightshift officer that the instrument didn't work when he attempted to use it.

When I arrived to look at the instrument it appeared to be working fine. I did a "test" and it came back with an Internal Standard Error message. I contacted the crime lab for further advice.

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: 6/29/20 ~ 8:35 am BMB 7/7/20

Out of Service Date/Time: 6/29/20 at about 1000 hours BMB 7/7/20 see info in status change about and additional notes below

For Use by SCDL

Additional Notes

I logged into the Juneau PD instrument remotely and checked the voltages and temperatures. Everything appeared normal. I ran a remote VOC and both checks of the internal standard passed. I contacted Sterling and told him we would monitor the instrument but that he should notify us if there were any further problems. The passing VOC put the instrument back in service.
BMB 7/7/20

Email completed form to dps.scdl.toxicology@alaska.gov

For questions contact the Breath Alcohol Program at 907-269-5740