

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Effective: 10/10/2022

Version: 3.0

Instrument #: 100687

Received From: AST/Galena

Date: 12/31/22

Analyst: C. O'Bryant

Reason for instrument return:

Removed from service for interference detected,
internal standard error + blank error.

Change instrument and case location in JusticeTrax Breath Alcohol Database

Ensure change in instrument status form is in JusticeTrax

Upload

Update agency information and polling status in DMHost Calling Program

Check for VOCs

From Intoximeters

Intoximeters Repair Paperwork (attached)

Screen was cracked upon return.
COB 12/31/22