

# Alaska Scientific Crime Detection Laboratory

## Change in Instrument Status Form

Issued: 6/11/2014  
Effective: 6/11/2014

Version: CIISF 2014 R0  
Status: Active

Date: 11/17/2016 Instrument Serial Number: 100364  
Instrument Location: ANCHORAGE JAIL  
Supervisor Name and Agency: Ofc P. GILBERT DSN 61033/APD

Describe the Status Change or Issue with the Instrument:

~~AN~~ OFFICER CONDUCTED A BREATH TEST ON A SUBJECT AND AFTER ALL CHECKS WERE DONE A MESSAGE POPULATED "ERROR UNDEFINED". I LOOKED ~~AT~~ AT THE LAST TEST DONE AND THE OFFICER'S TEST WAS NOT IN THE SYSTEM. I CONDUCTED A VOC AND AFTER ALL CHECKS WERE DOWN IT AGAIN POPULATED "ERROR UNDEFINED".  
I PLACED THE DMT "OUT OF SERVICE."

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: 11/18/16 00:20  
Out of Service Date/Time: 11/17/2016 (03:07:24)

For Use by SCDL

Additional Notes

Instrument restarted by Officer Gilbert followed by a non-drinking subject/test and VOC.  
Unable to recreate software message.  
Instrument returned to service  
OGB 11/18/16

Email completed form to [dps.scdl.toxicology@alaska.gov](mailto:dps.scdl.toxicology@alaska.gov)

For questions contact the Breath Alcohol Program at 907-269-5740

## **O'Bryant, Colleen S (DPS)**

---

**From:** Gilbert, Patrick E. <PGilbert@muni.org>  
**Sent:** Friday, November 18, 2016 2:19 AM  
**To:** O'Bryant, Colleen S (DPS)  
**Subject:** RE: 100364

Looks like the machine is good to go. I did the breath sample and all is working fine. A VOC was done as well.

**Officer Patrick Gilbert / DSN 61033**  
**Impaired Driving Enforcement Unit (IDEU)**  
**Days off MTW / (907)-729-7646**

---

**From:** O'Bryant, Colleen S (DPS) [<mailto:colleen.obryant@alaska.gov>]  
**Sent:** Thursday, November 17, 2016 7:19 AM  
**To:** Gilbert, Patrick E.  
**Cc:** Barnett, Brandi M (DPS); Foster, Charles R (DPS)  
**Subject:** 100364

**Patrick,**

**Thanks for the change in status form for 100364. I'd like to try some troubleshooting before we physically remove the instrument from service. Could you restart the instrument (turn off, count to 10, turn on) and attempt to perform a non-drinking subject test? Anyone can do this, doesn't have to be a breath test supervisor.**

**Colleen O'Bryant – Forensic Scientist III**  
**Blood and Breath Alcohol**  
**Alaska Scientific Crime Detection Laboratory**  
**4805 Dr. Martin Luther King Jr. Ave**  
**Anchorage, AK 99507**

**907-269-5783 (direct)**  
**907-350-2111 (cell)**

SUBJECT TEST

Alaska Department of Public Safety

DATAMASTER dmt: 100364

Date: 11/18/2016

Time: 01:53:30

SUBJECT NAME:

TEST TEST

LICENSE #: TEST

OPERATOR NAME:

PATRICK E GILBERT

OPERATOR'S #: 7367

DEPT/AGENCY:

ANCI

CASE/REPORT: 00-000000

TEST TYPE: TEST

EXTERNAL STANDARD INFORMATION

NOMINAL: 0.080

TARGET AT 30.04: 0.080

LOT #: AG608504

EXPIRATION: 03/25/2018

DIAGNOSTIC CHECK

BLANK TEST	PASSED	01:54
INTERNAL STANDARD	0.000	01:54
EXTERNAL STANDARD	VERIFIED	01:54
BLANK TEST	0.081	01:54
INTERNAL STANDARD	0.000	01:55
EXTERNAL STANDARD	0.000	01:56
BLANK TEST	0.000	01:57
INTERNAL STANDARD	0.082	01:57
EXTERNAL STANDARD	0.000	01:58

