

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Issued: 7/10/2015
Effective: 7/10/2015

Version: DMCIF 2015 R1
Status: Active

Instrument S/N 100363

Received From NPAS

Date 9/28/15

Analyst Colleen OByant

1. Reason for instrument return:

Instrument returned to NPAS for software issues.

2. Change instrument and case location in JusticeTrax Breath Alcohol Database

3. Ensure change in instrument status form is in JusticeTrax N/A

4. Upload N/A

5. Update agency information and polling status in DM Host Calling Program N/A

6. Check for VOCs N/A

From NPAS

7. NPAS Repair Paperwork (attached)

Work Ord No 49860

National Patent Analytical Systems, Inc.
Service Action Report

Serial No 100363

Date 8/18/15

Customer Alaska Dept of Public

Phone 907-269-5582

Contact Bolz

War Exp 1/13/13

Mx Exp

Ship Damage? Y N

Cust Pick Up?? Y N

Initial Comments: note any shipping damage or estimates required, along with accessories

SEE ATTACHED

Describe problem seen and/or attach documents received

multiple software issues

Keyword

Finding

Error Code

Enter Hours 1.0

Repair Comments and work Performed

Correctly updated software to v3.01 fixed sticking 5 way, Gave Int Std error. Recalibrated, ran standard tests ok

Solution Lot No 15040

| Voltage Set Values In | | | | Voltage Set Values Out | | | |
|-----------------------|--------|------------|-------|------------------------|--------|------------|-------|
| Lamp | 1.740 | Brth Vol | 1.580 | Lamp | 1.740 | Brth Vol | 1.540 |
| Meter | -0.221 | Filter Dif | 1.001 | Meter | -0.213 | Filter Dif | 1.014 |
| Stability | 0.002 | | | Stability | 0.001 | | |

| Part No | Qty Ord | Pick | Back Ordered | Total shin. | Description | Price Ea |
|---------|---------|------|--------------|-------------|---|----------|
| 9DMT | 1.0 | | | | DATAMASTER DMT SERVICE RETURN 100363 | 0.000 |
| NOTE | 1.0 | | | | NANUK CASE | |
| r500-1 | 1.0 | | | | Service Charge for Repair | 119.000 |

Date Completed 9/1/15

Technician r. martin

Pick up Signature