

Alaska Scientific Crime Detection Laboratory

Change in Instrument Status Form

Issued: 6/11/2014
Effective: 6/11/2014

Version: CIISF 2014 R0
Status: Active

Date: 10/14/16 Instrument Serial Number: 100364

Instrument Location: Anchorage Jail

Supervisor Name and Agency: Brandi Barnett

Describe the Status Change or Issue with the Instrument:

Instrument was giving blank error + internal standard status messages.

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: _____

Out of Service Date/Time: _____

For Use by SCDL

Additional Notes

See attached email. I ran a VOC successfully after purging the sample chamber on 10/17/16. It appeared the problem was environmental. I requested that Officer Palmer verify the vent tube was still attached which he confirmed on 10/14/16. 10/13/16. BMB. 10/19/16
BMB 10/14/16

Email completed form to dps.scdl.toxicology@alaska.gov

For questions contact the Breath Alcohol Program at 907-269-5740

Barnett, Brandi M (DPS)

From: Palmatier, Kevin M. <KPalmatier@muni.org>
Sent: Thursday, October 13, 2016 10:58 PM
To: Barnett, Brandi M (DPS)
Subject: Re: Datamaster 100364
Attachments: image1.JPG

Everything appears to be hooked up fine.
[cid:DA610117-B08B-44A5-8BB1-C4B6DBAF74E4]
Kevin Palmatier

Sent from my iPhone

On Oct 11, 2016, at 10:29 AM, Barnett, Brandi M (DPS)
<brandi.barnett@alaska.gov<<mailto:brandi.barnett@alaska.gov>>> wrote:

Did you get a chance to look at the vent tube yet?

Brandi

From: DPS SCDL Toxicology (DPS sponsored)
Sent: Friday, October 07, 2016 7:00 AM
To: Foster, Charles R (DPS); Barnett, Brandi M (DPS); O'Bryant, Colleen S (DPS)
Subject: FW: Datamaster 100364

From: Palmatier, Kevin M. [KPalmatier@muni.org<<mailto:KPalmatier@muni.org>>]
Sent: Friday, October 07, 2016 5:24 AM
To: DPS SCDL Toxicology (DPS sponsored)
Subject: Datamaster 100364

Good morning,

I just wanted to let you know that I had some issues with Datamaster 100364 that is at the Anchorage Jail today. I was trying to get a breath sample from a subject and after capturing the breath sample, the Datamaster displayed a Blank Error. I moved the subject away from the Datamaster and set it up for another sample. I received this error two more times. I was able to find a fan and attempted to clear out the Datamaster room. I set it up again and received an internal standard error. I changed rooms and tested the subject on another Datamaster without any issues.

Please let me know if there is anything I need to do with the above listed Datamaster.

Thanks,
Ofc Palmatier

MUNICIPALITY OF ANCHORAGE
Police Department

OFFICER PALMATIER
Impaired Driving Enforcement Unit

4501 Elmore Road
Anchorage, Alaska 99507-1599

Telephone (907) 786-8771
Dispatch (907) 786-8900

Email: kpalmatier@muni.org<mailto:kpalmatier@muni.org>

