

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Effective: 11/10/2021

Version: 2.0

Instrument #: 100384

Received From: Sand Point PD

Date: 8/16/22

Analyst: C. O'Bryant

Reason for instrument return:

DMT Returned to reconcile data. DMT arrived with external standard still installed and missing power cord + key board.

Change instrument and case location in JusticeTrax Breath Alcohol Database

Ensure change in instrument status form is in JusticeTrax

Upload

Update agency information and polling status in DMHost Calling Program

Check for VOCs

From Intoximeters

Intoximeters Repair Paperwork (attached)

* Upon start up the instrument displayed a white screen. I opened the cover + unplugged + plugged back in the J24 cable connected to the screen + controller board. The screen started up + functioned normally following this.