VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Serial #: 100695 Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program Date:01/01/2013

External Standard Test Values

EXTERNAL STANDARD INFORMATION

NOMINAL: 0.080

TARGET AT 29.51: 0.078 LOT #: 22612080A1 EXPIRATION: 09/01/2014

BLANK TEST	0.000	02:48
INTERNAL STANDARD	VERIFIED	02:48
EXTERNAL STANDARD	0.079	02:48
BLANK TEST	0.000	02:49
EXTERNAL STANDARD	0.079	02:50
BLANK TEST	0.000	02:50
EXTERNAL STANDARD	0.079	02:51
BLANK TEST	0.000	02:52
EXTERNAL STANDARD	0.079	02:52
BLANK TEST	0.000	02:53
EXTERNAL STANDARD	0.078	02:53
BLANK TEST	0.000	02.54

Average = 0.0788Std Dev = 0.0004

Diagnostic Check

VERSTONS DMT: 2.00 PIC: 3.03 Modem: 2.3 Questions: 2.2

TEMPERATURES Sample Chamber = 48.5°C PASSED Breath Tube = 41.6°C PASSED PUMP INFO Flow Rate = 3.976 L/MPASSED DETECTOR INFO PIIMP ON PASSED PUMP OFF PASSED FILTER INFO Filter 1 PASSED Filter 2 PASSED Filter 3 PASSED INTERNAL STANDARD PASSED

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Scientific Director State Breath Alcohol Program

Subscribed and sworn before me this 31 day of

Nikki Roth, Notary Public My Commission Expires With Office



