

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Issued: 7/10/2015
Effective: 7/10/2015

Version: DMCIF 2015 R1
Status: Active

Instrument S/N 100411

Received From NPAS

Date 7/30/15

Analyst Brandi Barnett

1. Reason for instrument return:

Return from NPAS

2. Change instrument and case location in JusticeTrax Breath Alcohol Database

3. Ensure change in instrument status form is in JusticeTrax N/A

4. Upload N/A

5. Update agency information and polling status in DM Host Calling Program

6. Check for VOCs N/A

From NPAS

7. NPAS Repair Paperwork (attached)

Packing List
From: National Patent Analytical Systems Inc.
 2090 Harrington Mem. Rd., Mansfield, Ohio 44903

Cust No. AK6056A **Ord Date** 6/18/15
Sales Order No. 49663

Ship Date 7/14/15
Packing List No. 74288

Cust Purch Ord No. COLLEEN

Ship To:
 Alaska Dept of Public Safety
 Scientific Crime Detection
 4805 Dr. Martin Luther King Dr.

Bill To:
 Alaska Dept of Public Safety
 Scientific Crime Detection
 4805 Dr Martin Luther King Dr.

Anchorage Ak 99507

Anchorage Ak 99507 **Freight**

25.00

Attn: Ship Via FedEx

Comments Service for100411

Part / Drawing No / Rev No	Quantity Ordered	This Shipment	Back Ordered	Total Ship To Date	Description	NP Invent Ref
No Draw # / NONE 100411 Serial Nos.	1.0	1		1.0	DATAMASTER DMT SERVICE RETURN	9DMT
No Draw # / NONE Serial Nos.	1.0	1		1.0	NANUK CASE	NOTE
No Draw # / NONE Serial Nos.	1.0	1		1.0	Service Charge for Repair	r500-1
No Draw # / NONE Serial Nos.	1.0	1		1.0	Service Charge for Repair, Additional Hours	r500-2

Work Ord No 49663

National Patent Analytical Systems, Inc.
Service Action Report

Serial No 100411

Date 6/18/15

Customer Alaska Dept of Public

Phone 907-269-5582

Contact Bolz

War Exp 5/4/13

Mx Exp

Ship Damage? Y N

Cust Pick Up?? Y N

Initial Comments: note any shipping damage or estimates required, along with accessories

Describe problem seen and/or attach documents received i/f, int std error, fatal errors

Keyword

Finding

Error Code

Enter Hours 2.0

Repair Comments and work Performed

Updated software to address errors, inspected unit, set voltages, calibrated, ran all standard tests ok

Solution Lot No 15040

Voltage Set Values In		Voltage Set Values Out	
Lamp	1.770	Brth Vol	1.570
Meter	0.360	Filter Dif	0.883
Stability	0.002	Lamp	1.840
		Brth Vol	1.560
		Meter	-0.024
		Filter Dif	0.958
		Stability	0.002

Part No	Qty Ord	Pick	Back Ordered	Total sbio	Description	Price Ea
9DMT	1.0				DATAMASTER DMT SERVICE RETURN 100411	
NOTE	1.0				NANUK CASE	
r500-1	1.0				Service Charge for Repair	0.000
r500-2	1.0				Service Charge for Repair, Additional Hours	0.000

Date Completed 7/8/15

Technician R. MARTIN

Pick up Signature