

# Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form

Issued: 7/10/2015  
Effective: 7/10/2015

Version: DMCIF 2015 R1  
Status: Active

Instrument S/N 100852

Received From NPAS

Date 9/29/15

Analyst Colleen O'Bryen

1. Reason for instrument return:

Returned from NPAS for software issues

2. Change instrument and case location in JusticeTrax Breath Alcohol Database

3. Ensure change in instrument status form is in JusticeTrax  N/A

4. Upload  N/A

5. Update agency information and polling status in DM Host Calling Program  N/A

6. Check for VOCs  N/A

From NPAS

7. NPAS Repair Paperwork (attached)

Work Ord No 49916

National Patent Analytical Systems, Inc.  
Service Action Report

Serial No 100352

Date 9/2/15

Customer Alaska Dept of Public

Phone 907-269-5582

Contact Bolz

War Exp 1/13/13

Mx Exp

Ship Damage?  Y  N

Cust Pick Up??  Y  N

Initial Comments: note any shipping damage or estimates required, along with accessories

Describe problem seen and/or attach documents received

FATAL ERRORS

Keyword

Finding

Error Code

Enter Hours 1.0

Repair Comments and work Performed

Ticket DB error. Couldnt update from USB, wiped software, cleared registry, loaded software, calibrated, ran standard tests ok

Solution Lot No 15040

Voltage Set Values In				Voltage Set Values Out			
Lamp	1.720	Brth Vol	1.540	Lamp	1.800	Brth Vol	1.550
Meter	0.340	Filter Dif	0.789	Meter	-0.193	Filter Dif	0.881
Stability	0.002			Stability	0.002		

Part No	Qty Ord	Pick	Back Ordered	Total shin	Description	Price Ea
9DMT	1.0				DATAMASTER DMT SERVICE RETURN 100352	
NOTE	1.0				NANUK CASE	
R500-1	1.0				Service Charge for Repair	119.000

Date Completed 9/18/15

Technician R. MARTIN

Pick up Signature