Alaska Scientific Crime Detection Laboratory
DataMaster Evaluation Form
Effective: 10/10/2022 Version: 3.0
Instrument #: $100682$ Analyst(s): $COBrant$ Start Date: $N/6/22$ Date Received: $N/6/22$
Reason for instrument evaluation Unable to load Hicket database
Visual Inspection Clean Clean Replace missing screws
Cermetex Box (modem) is secure
Double Check Valve secured to wall
O-ring inspection and ensure 2 replacement O-rings in bag
As Found-Conditions: Barometric Pressure (inHg):Temperature (C):Relative Humidity (%) Linearity Test (attached) 0.020g/210L Mnf:Lot:Exp: 0.080g/210L Mnf:Lot:Exp: 0.300g/210L Mnf:Lot:Exp: Notes Notes
<ul> <li>VOC (attached)</li> <li>Verify software version</li> <li>VOC Internal Standard value ≤ 3.00%</li> <li>Non-Drinking Subject Test (attached)</li> </ul>
Condition of item: Good Meeds Repair
Ready for Adjustment     Ready for Calibration
End Date:
Page 1 of 1
All printed copies are uncontrolled. Approved by: Chemistry Supervisor