

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Issued: 2/2/2015
Effective: 2/2/2015

Version: DMCIF 2015 R0
Status: Active

Instrument S/N 100346
1006-03042515

Received From AST/Soldotna

Date 6/25/15

Analyst Colleen OBryant

1. Reason for instrument return:

USB connectivity issues as well as uploading
issues, gave message "external storage device not
found"

2. Change instrument and case location in JusticeTrax Breath Alcohol Database

3. Ensure change in instrument status form is in JusticeTrax

4. Visual Inspection

5. Upload

6. Update agency information and polling status in DM Host Calling Program

7. Check for VOCs

From NPAS

8. NPAS Repair Paperwork (attached)