

Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form

Effective: 11/10/2021

Version: 2.0

Instrument #: 100414

Analyst(s): Carlie Glaister

Start Date: 2/24/2022

Reason for instrument evaluation

Screen is inoperable

Visual Inspection

Clean

Replace missing screws

Cermetex Box (modem) is secure

Ensure modem cable present

Double Check Valve secured to wall

All tubing and connections secure

O-ring inspection and ensure 2 replacement O-rings in bag

~~connections to not secured.~~

*OK 2/24/22
screen*

Troubleshooting/Repair

Secured double check valve to wall with double sided tape.
Removed screen.

Screen to be sent to factory for repair. No screens available at this time for replacement.

Add replacement screen when they arrive.

VOC (attached)

Verify software version

VOC Internal Standard value $\leq 3.00\%$

Non-Drinking Subject Test (attached)

OK 2/24/2022

Return to Factory

Ready for Calibration

Ready for Certification

Proceed to Close Out

Further Evaluation Needed

End Date: 2/24/2022