## Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form Fifective: 10/10/2022 Version: 3.0 Instrument #: 100377 Received From: Inteximeters Analyst: Kaija Svisher Supervised casework with Derek Walton 06/14/2024 Reason for instrument return:

Repair. Filter wheel errors. Problem with uploading.	
Change instrument and case location in JusticeTrax Breath Alcohol Database	
insure change in instrument status form is in JusticeTrax	
Instrument inadvertently uploaded after reffrom Intox. This has no impact to performa DJW 06/14/2024	turn ince
Ipdate agency information and polling status in DMHost Calling Program 🗹	
Check for VOCs 🔽	

From Intoximeters

Intoximeters Repair Paperwork (attached)



Experience · Service · Integrity

## Service Record

SRO Number: SRO-060460

Received Date: 01/17/2024

Completion Date: 05/07/2024

LEGACY DMT

SRO Type: REPAIR

SRO Description: DMT Repair

**Customer ID:** C000003651

Customer PO: 01172024-OBryant

Contact: Colleen OBryant

Email: Colleen.obryant@alaska.gov

LEGDMT

Ship To:

Units on SRO

Alaska Scientific Crime Detection Laboratory 4805 Dr MLK Jr Ave

Anchorage AK 99507 USA

Service Line:	**All Instruments are Calibrated to Factory Specifications**						
Unit Type:	(LEGDMT) LEGACY DMT						
Serial Number:	50100377						
Reason:	Reason: DMT Status Messages						
	FitWl er						
Reason Notes:	Customer Reported Issue: Filter wheel errors, won't upload						
Resolution:							
	RPL comp						
Resolution Notes:	tes: Replace filter wheel solenoid.						
Parts:	Qty	Part Number	Part Description				
	1	59-0080-00	DMT ELEC ASSY SOLENOID FILTR WHEEL 41344				
Reason:	DMT Functional						
	No Displ						
Reason Notes:	Touch screen random not operational.						
Resolution:							
	RPL comp						
Resolution Notes:	Replace the touch screen membrane.						
	1	001199-00	DMT ELEC REGULATOR V LDO 3.3V OUT 30V IN				
Reason:	Upgrade to newer version						
	Upgrade						
Reason Notes:	Digital pot not complete.						
Resolution:	: Update						
	DMTVR3						
Resolution Notes	Digital pot mode is complete. (No charge for this)						
	h .	59-0101-00	DMT ELEC MEMBRANE 8.4IN TOUCH REPL 37553				

## **Order Packing Slip Report**

Intoximeters
Experience · Service · Integrity

Pack Date: 5/21/2024

5/21/2024 7:30:30 AM

Packing Slip: 307654

SRO: SRO-060460

): SRO-000400

From: Intoximeters, Inc

2081 Craig Road St Louis MO 63146

USA

**Customer:** C000003651 Ship To: 1

Alaska Department of Public Safety

4805 Dr MLK Jr Ave Anchorage AK 99507

**USA** 

Alaska Scientific Crime Detection

Laboratory

4805 Dr MLK Jr Ave Anchorage AK 99507

USA

Customer PO: CC-01172024-

**OBryant** 

Ship Via: FDX

**FED EX** 

**Weight: 36.00** 

Delivery Terms: 2D

Shipping Terms: ORIGIN

# Packages: 1

**Customer Ship Account:** 

Line Number	Item	UM	Ordered Quantity	Shipment Quantity	Backorder Quantity					
1	LEGDMT	LEGACY	DMT		50100377					
1	59-0080-00	EA		1.00	0.00					
DMT ELEC ASSY SOLENOID FILTR WHEEL 41344										
2	001199-00	EA		1.00	0.00					
DMT ELEC REGULATOR V LDO 3.3V OUT 30V IN										
3	59-0101-00	EA		1.00	0.00					
DMT ELEC N	MEMBRANE 8	.4IN TOUC	H REPL 37553							