

Alaska Scientific Crime Detection Laboratory

Change in Instrument Status Form

Issued: 6/11/2014
Effective: 6/11/2014

Version: CIISF 2014 R0
Status: Active

Date: 3/2/16 Instrument Serial Number: 100695

Instrument Location: Wasilla PD

Supervisor Name and Agency: Charles Foster - SCDL

Describe the Status Change or Issue with the Instrument:

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: _____

Out of Service Date/Time: _____

For Use by SCDL

Additional Notes

Email documenting troubleshooting performed by Colleen O'Bryant on 2/18/16 is attached.

Email completed form to dps.scdl.toxicology@alaska.gov

For questions contact the Breath Alcohol Program at 907-269-5740

From: O'Bryant, Colleen S (DPS)
Sent: Monday, February 29, 2016 7:19 AM
To: Foster, Charles R (DPS)
Subject: Fwd: datamaster recert

On 2/18/16 I went to WPD after reported USB connectivity issues. I restarted the instrument and the printer. Unplugged all exterior USB connections and reconnected. I didn't not experience and screen or keyboard issues.

I initiated a VOC and left the instrument.

Colleen

Sent from my iPhone

Begin forwarded message:

From: "O'Bryant, Colleen S (DPS)" <colleen.obryant@alaska.gov>
Date: February 17, 2016 at 11:37:19 AM AKST
To: Michael Bonadurer <mbonadurer@ci.wasilla.ak.us>
Subject: Re: datamaster recert

The USB pop up is not uncommon on start up. Remember to press the X, not the OK.
Locks up when pressing run? Is the screen being unresponsive? Does a restart fix it? What frequency is this occurring?

Colleen

Sent from my iPhone

On Feb 17, 2016, at 9:27 AM, Michael Bonadurer <mbonadurer@ci.wasilla.ak.us> wrote:

Hey Colleen,

Ben Dudley informed me we are having issues with our datamaster at WPD. It locks up when you hit the run button and when you restart it a USB error window pops. Thoughts, feelings?

Mike