

Alaska Scientific Crime Detection Laboratory

Change in Instrument Status Form

Issued: 6/11/2014
Effective: 6/11/2014

Version: CIISF 2014 R0
Status: Active

Date: 9/9/15 Instrument Serial Number: 100359

Instrument Location: Dillingham PD

Supervisor Name and Agency: Dan Pasquariello (PPD) / C.O. Bryant (SCDL)

Describe the Status Change or Issue with the Instrument:

Chief Pasquariello contacted the lab regarding a 'fatal system error' message an operator got while trying to set up a subject test. When using the drop down menu to select the test type, the DMT would give a 'fatal system error' and restart. This occurred multiple times. The operator attempted to set up the subject test again without using the drop down menu and did not receive the message. The message was not occurring otherwise. COB 9/9/15

If instrument was placed in service or removed from service, fill out the appropriate sections below.

In Service Date/Time: _____

Out of Service Date/Time: _____

For Use by SCDL

Additional Notes

I reinstalled 3.01 software via the modem line, and performed a VOC remotely. I requested the chief perform tests attempting to recreate the (problem) and he was unable to do so. The instrument was left in service. (See attached email) COB 9/9/15

Email completed form to dps.scdl.toxicology@alaska.gov

For questions contact the Breath Alcohol Program at 907-269-5740

O'Bryant, Colleen S (DPS)

From: Daniel Pasquariello <chiefofpolice@dillinghamak.us>
Sent: Tuesday, September 08, 2015 6:01 PM
To: O'Bryant, Colleen S (DPS)
Subject: RE: Installed 3.01 software again

Colleen,

No errors. Thank you.

Chief Daniel Pasquariello
Dillingham Dept. of Public Safety
PO Box 869
Dillingham, Alaska 99576
(907) 842-5354
(907) 842-5785 fax
chiefofpolice@dillinghamak.us



From: O'Bryant, Colleen S (DPS) [mailto:colleen.obryant@alaska.gov]
Sent: Tuesday, September 08, 2015 1:33 PM
To: chiefofpolice@dillinghamak.us
Subject: Installed 3.01 software again

Please attempt to recreate the problem by running a few tests and let me know if you get a fatal system error.

Thanks,

**Colleen O'Bryant – Forensic Scientist III
Blood and Breath Alcohol
Alaska Scientific Crime Detection Laboratory
4805 Dr. Martin Luther King Jr. Ave
Anchorage, AK 99507**

**907-269-5783 (direct)
907-350-2111 (cell)**