| Alaska Scientific Crime Detection Laborato | ry |
|--|----|
|--|----|

DataMaster Evaluation Form

Issued: 7/10/2015 Effective: 7/10/2015 Version: DMEF2015 R0 Status: Active

| Instrum | nent #_ | 100689 | |
|---------|---------|--------|---|
| Date: _ | 101 | 23/15 | _ |

Analyst: Brandi Barnett

- 1. Reason for instrument evaluation Proventative Maintenance
- 2. Visual Inspection
- 3. Troubleshooting/Repair

- 4. VOC (attached)
- Internal Standard value ≤ 3.00% Yes No
 Verify software version
- Return to Factory for Repair 🔲
 - Ready For Certification

Ready For Calibration

VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument

State of Alaska Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Serial #: 100689

Date:10/23/2015

| External Standard Test Values | Diagnostic Check |
|--|--|
| EXTERNAL STANDARD INFORMATION NOMINAL: 0.080 TARGET AT 29.70: 0.079 LOT #: 01315080A2 EXPIRATION: 02/05/2017 TANK PRESSURE: 237 psi BLANK TEST 0.000 13:46 | VERSIONS DMT 3.01 PIC: 93 Modem: 2.6 Questions: 2.2 TEMPERATURES |
| INTERNAL STANDARD VERIFIED 13:46 Xq = 0.106 (0.67%) EXTERNAL STANDARD 0.678 13:47 | Sample Chamber = 48.9° C PASSED (44.0°C - 52.0°C) |
| X[1] = 0.0784 (0.0000) (0.0020) BLANK TEST 0.000 13:47 EXTERNAL STANDARD 0.078 13:48 X[1] = 0.0782 (-0.0006) (0.0020) | Breath Tube = 48.1°C PASSED (38.0°C - 50.0°C) SETTINGS |
| BLANK TEST 0.000 13:49 EXTERNAL STANDARD 0.078 13:49 | Lamp Voltage = 1.67 V PASSED (1.05 V - 2.45 V) |
| $X[1] = 0.0784 \{0.0001\} \{0.0015\}$ BLANK TEST 0.000 13:50 | Cooler Voltage = 1.75 V PASSED (0.99 V - 2.31 V) |
| EXTERNAL STANDARD 0.078 13:50 X[1] = 0.0782 (0.0001) (0.0010) | $\begin{array}{rcl} \text{Bias Voltage} &= 80 \text{ V} & \text{PASSED} \\ (49 \text{ V} - 115 \text{ V}) \end{array}$ |
| BLANK TEST 0.000 13:51 EXTERNAL STANDARD 0.078 13:51 X[1] = 0.0782 (0.0002) (0.0022) BLANK TEST 0.000 13:52 | Chopper Freq = 555 Hz PASSED (475 Hz - 575 Hz) Barometer = 29.7 in |
| Average = 0.0780 Std Dev = 0.0000 | PUMP INFO Flow Rate = 4.603 L/M PASSED (3.500 L/M - 6.500 L/M) |
| | DETECTOR INFO PUMP ON PASSED (0.002 V <= 0.010 V) |
| | PUMP OFF PASSED (0.002 V <= 0.010 V) |
| | FILTER INFO Filter 1 PASSED |
| | Filter 2 PASSED |
| | Filter 3 PASSED |
| | $\begin{array}{rcl} & & & & & \\ & & & & & \\ & & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & & & \\ & & & & & \\ & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & &$ |
| | Ag = 0.100 0.425 0.000 4.000 |

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Nita J.Bolz Scientific Director State Breath Alcohol Program

Subscribed and sworn before me this _____ day of _____, 20 _____