## VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Date:06/10/2016

## External Standard Test Values

EXTERNAL STANDARD INFORMATION

NOMINAL: 0.080

Serial #: 100698

TARGET AT 29.78: 0.080 LOT #: 01315080A2 EXPIRATION: 02/05/2017 TANK PRESSURE: 1011 psi

0.000	22:30
VERIFIED	22:30
0.080	22:30
0.000	22:31
0.081	22:31
0.000	22:32
0.081	22:32
0.000	22:33
0.080	22:33
0.000	22:34
0.080	22:34
0.000	22:35
	VERIFIED 0.080 0.000 0.081 0.000 0.081 0.000 0.080 0.080 0.000

Average = 0.0804Std Dev = 0.0005

## Diagnostic Check

VERSIONS DMT: 3.02 PIC: 3.03 Modem: 2.6 Questions: 2.2

TEMPERATURES	
Sample Chamber = 48.9°C	PASSED
Breath Tube = 46.6°C	PASSED
PUMP INFO	
Flow Rate = 4.589 L/M	PASSED
DETECTOR INFO	
PUMP ON	PASSED
PUMP OFF	PASSED
FILTER INFO	
Filter 1	PASSED
Filter 2	PASSED
Filter 3	PASSED
INTERNAL STANDARD	PASSED

I, Charles R. Foster, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

7/25/16

Charles R. Foster Scientific Director

State Breath Alcohol Program

Subscribed and sworn before me this aday of July

Nikki Roth, Notary Public My Commission Expires With Office





8/1/22/10